S. PneumoStrip

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S. PneumoStrip

Test for the identification of 76 serotypes of Streptococcus pneumoniae

INTENDED USE

The S. PneumoStrip test is a test based on the reverse blot technique that allows the identification, in base of its genetic sequence, of 76 serotypes of *S. pneumoniae*:

- Strip A: 1, 3, 4, 5, 6A, 6B, 6C, 6D, 7F/7A, 9A/9V, 14, 18A, 18B/18C, 18F, 19A, 19F, 23A, 23B, 23 F.
- Strip B: 2, 8, 9N/9L, 10A, 10B, 10F/C, 11A/D, 11B, 11C, 11F, 12A/46, 12B/44, 12F, 15B/15C, 17F, 20, 22F/22A, 33F/33A, 37.
- Strip C: 7B, 7C/40, 15A, 15F, 16F, 19B/19C, 21, 25A/25F, 38, 24A, 24B/24F, 31, 32A/32F, 33B/33D, 33C, 35A, 35C, 35F, 47F, 41A, 41F.

Of the 76 serotypes, 42 are identified individually and 34 as pairs (for example, the presence of serotypes 7F or 7A).

Among the known serotypes are the 23 included in the various vaccines available (7-valent, 10-valent, 13-valent and 23-valent):

Vaccine	Serotypes <i>S. pneumoniae</i>
PCV7	4, 6B, 9V, 14, 18C, 19F, 23F
PCV10	1, 4, 5, 6B, 7F, 9V, 14, 18C, 19F, 23F
PCV13	1, 3, 4, 5, 6A, 6B, 7F, 9V, 14, 18C, 19A, 19F, 23F
PCV23	1, 2, 3, 4, 5, 6B, 7F, 8, 9N, 9V, 10A, 11A, 12F, 14, 15B, 17F, 18C, 19A, 19F, 20, 22F, 23F, 33F.

PRINCIPLE

The S. PneumoStrip kit is based on the principle of reverse hybridization, and allows the detection and identification of DNA samples, extracted from bacterial culture, of 42 serotypes of *S. pneumoniae*.

Infections of *S. pneumoniae* represent an important public health problem in both industrialized nations as well as those that are less-developed. It is responsible for high mortality and lethality, as it is one of the main causative agents for a large variety of clinical conditions, from benign infections like otitis media and acute sinusitis to more severe infections such as septicaemia, meningitis and pneumonia.

Among the interesting structures of this bacteria is a capsule outside the cell wall, made of complex polysaccharides. This capsule protects the bacterium from phagocytosis and is essential to its virulence. In fact, the non-encapsulated forms of the bacteria are, normally, not virulent.

The antigenic composition of the capsule is variable between the various strains and allows *S. pneumoniae* to be grouped into over 100 different capsular serotypes and approximately 45 serogroups. Serotypes that present cross-immunogenicity are considered to belong to the same serogroups, for example, 6A and 6B. The identification of each serotype is made using an antigen-

antibody reaction with specific antisera, which results in swelling of the capsule, a phenomenon known as “Quellung”, a German term that means swelling.

The procedure of *S. PneumoStrip* kit consists in three steps: a) DNA extraction, b) PCR amplification and c) Hybridization/developing.

a) DNA extraction

Reagents not included in the kit.

Any standard DNA purification protocol can be used, including rapid DNA extraction from colonies by boiling (see “SAMPLES” section).

b) PCR amplification

The *S. PneumoStrip* kit includes the reagents needed to perform specific genetic sequence amplification of the serotypes detected with the kit.

For the differentiation of the various serotypes of *S. pneumoniae* by PCR, the *cps* locus is used, responsible for the synthesis of the bacterial capsid. This locus shows a common genetic organisation across all serotypes: flanked by two conserved genes, *dex B* and *aliA*, it includes four primary regulatory and processing genes (*wzg*, *wzh*, *wzd* and *wze*), highly conserved in nearly all serotypes, and a central portion that contains serotype-specific genes (such as *wzy* and *wzx*) that code for proteins such as polysaccharide synthetase, flippase, various glycosyl transferases, etc. These are specific genes that are used to distinguish one serotype from another.

During PCR, three control genes are also amplified:


- Autolysin, encoded by the gene *lytA*, a protein necessary for the pathogenesis of *S. pneumoniae* is a well-characterised, highly-conserved virulence marker, that allows the bacteria to be distinguished from other genotypically similar species such as *S. mitis* and *S. orale* (McAvin J. et al., 2001).
- Pneumolysin, a pore-forming toxin 53 kDa long, encoded by the *ply* gene, is another key protein in *S. pneumoniae* virulence, present in all the described serotypes of the bacteria (Hirst R.A. et al., 2004).
- A gene of the *cps* locus, *CpSA*, present in all pneumococci with capsids. The expression of this capsid is important for the survival of the bacteria in blood and is strongly associated with the ability of the pneumococcus to develop an invasive infection (Bentley S.D. et al., 2006).

c) Hybridization and development

Specific binding of PCR-amplified DNA fragments takes place during this step, using a series of probes that are covalently bound to a nylon membrane.

The membranes are covalently bound to probes specific to each detected serotype, along with probes for the control genes (*CpSA*, *ply* and *lyt*), a visible control line probe and between one and three lines, black or red, for the control of the strip's position to aid in the interpretation of results (see attached drawing).

The detection of the fragments that hybridize to the different probes is carried out using a streptavidin-peroxidase conjugate that binds to a biotin marker added to the PCR-amplified DNA fragments. Following the addition of a peroxidase substrate (TMB), a blue precipitate appears in the locations where hybridization has taken place. The end result of the test is a pattern of bands that can be interpreted with the help of a control strip.

KIT CONTENTS

S. PneumoStrip			Kit 16 tests
	Membranes	STRIPS	16 + 16 + 16 tiras
	8 channel tray	PL	6
Reactivos de PCR REAG PCR	PCR Mix	REAG PCR MIX	0,8 ml
	Primers	REAG PRIMER	0,06 + 0,06 + 0,06 ml
	Taq	REAG TAQ HS	0,030 ml
	Denaturing	SOLN DN	1 ml
	Hybridization buffer	BUF HYB	125 ml
	Wash 1 buffer	BUF WASH 1	250 ml
	Conjugate	CONJ HRP	125 ml
	Wash buffer 2	BUF WASH 2	350 ml
	Substrate	SUBS TMB	65 ml
Instructions	Instructions		1
	Interpretation chart		1 + 1 + 1

MATERIALS NOT INCLUDED IN THE KIT

The following additional material is required when using the kit:

1. Microtubes for PCR
2. Micropipettes and micropipette tips (sterile or UV-irradiated and ideally with a filter)
3. Tweezers and a pencil (optional)
4. Chronometer
5. Thermometer

REQUIRED EQUIPMENT FOR KIT DEVELOPMENT

1. Thermocycler
The following thermocyclers have been successfully used with the kit: 2720 Thermal Cycler (Applied Biosystems), SimpliAmp Thermal cycler (Applied Biosystems), MJ Mini Gradient Thermal Cycler (BioRad), Mastercycler Personal (Eppendorf), S1000 Thermal Cycler (BioRad).
2. Thermoshakers
The kit has been validated with the following thermoshakers: PST-60HL (Biosan S.L.), Profiblot T30 AND T48 (Tecan), Autoblot 3000H (medTEC), TENDIGO (Fujirebio), Dynablot Heat (Dynex).
3. BLOTriX S1 or BLOTriX R2 scanner (BioSciTec, optional).

PRECAUTIONS

1. Only use the reagents *in vitro*.
2. Strictly follow the instructions provided. Modifying any of the prescribed **steps or temperatures** may severely affect the test results.
3. All the materials and reagents used must be free of DNAses. The use of filtered pipette tips is recommended for PCR preparation to prevent aerosol contamination. Ensure that all the usual precautions are taken in the preparation of the amplifications. Use autoclaved material for the hybridization/development process.
4. Store the kit components as indicated in the instructions.
5. Do not exchange components from kits with different lot numbers.
6. Do not use kit components after the expiration dates.
7. Strips are for single-use.
8. Supplied trays are for single-use.
9. Each channel of the supplied trays is for just one strip.
10. If the package is broken, the product can still be used providing none of the components have been damaged.
11. The used product should be discarded in compliance with current legislation.
12. Patient samples must always be treated as potentially infectious. Environmental and safety standards must be adhered to.
13. Once it is used, it is advisable to manipulate the strip taking into account the same considerations as for the sample. The strips have to be manipulated with the proper precautions and should be managed as biohazardous material.
14. The denaturation solution contains < 2% NaOH and is an irritant to both eyes and skin (H314 and P280, P305, P351, P338, P310).
15. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the CE marking and component lots.

STORAGE

Store all reagents at 2-8 °C. Precipitates may appear in some solutions (hybridization buffer, wash buffers 1 and 2, conjugate) owing to storage at 2-8 °C. The precipitates are dissolved again when they are brought back to room temperature (conjugate) or heated to 42 °C (hybridization buffer and wash buffers 1 and 2).

The expiry date of all the reagents is printed on the label.

SAMPLES

The test has been designed and validated for its use with DNA obtained from bacterial culture (for example, in Blood Agar plates).

Some of the DNA extraction methods that have been evaluated with satisfactory results with S. PneumoStrip kit are:

- DNA extraction with magnetic particles (NucliSENS easyMAG; BioMérieux). Dilution of DNA to 1 ng/ μ L.
- DNA extraction with silica resins (QIAamp DNA Mini Kit; Qiagen). Dilution of DNA to 1 ng/ μ L.
- Rapid DNA extraction in TE buffer: an isolated colony or equivalent of culture is introduced into a microtube with 50 μ L of PCR-quality TE buffer. It is heated to 100 °C for 15 minutes and centrifuged for 2 min at 14,000 rpm. 5 μ L of the supernatant are amplified.

Store the DNA samples at 2-8 °C if they are going to be used shortly after, or at -20 °C for longer periods of storage.

S. PNEUMOSTRIP PROCEDURE

1.- Polymerase chain reaction

PCR preparation

Important: before opening the vials with the PCR reagents centrifuge them briefly. This will ensure that all the contents will be at the bottom of the tube.

- Prepare the required PCR tubes based on the number of DNAs to be amplified.
- Add to each PCR tube: 36,5 μ L of PCR premix + 2,5 μ L of primers A + 2,5 μ L of primers B + 2,5 μ L of primers C + 1 μ L of Taq + 5 μ L of DNA (1 ng/ μ L). Mix well. If possible, keep the reagents and the mixtures at 2-8 °C during the preparation.
If various DNA samples are to be amplified, it is recommended to prepare a common mixture with all the reagents to finally add 45 μ L of mixture + 5 μ L of DNA. For example:

N° of PCRs	PCR premix	Primers A	Primers B	Primers C	Taq
1	36,5 μ L	2,5 μ L	2,5 μ L	2,5 μ L	1 μ L
3	146 μ L	10 μ L	10 μ L	10 μ L	4 μ L
5	219 μ L	15 μ L	15 μ L	15 μ L	6 μ L
8	328,5 μ L	22,5 μ L	22,5 μ L	22,5 μ L	9 μ L

** The mixtures containing all PCR reagents should always be prepared in excess to compensate for the loss of volume that takes place during the pipetting process.*

Due to requirements of "good laboratory practice", the user should also include a negative control (e.g. water or TE buffer) to exclude contamination and, if required, a positive control (not included in the kit).

Amplification

Regard the importance of getting a suitable concentration of the DNA to get a proper result.

Place the tubes into the thermocycler (if required, add 1 drop of Nujol oil on top of every PCR reaction) and amplify DNA by the following program:

96 °C, 5 min
40 cycles of: 96 °C, 15 sec
58 °C, 1 min
72 °C, 30 sec
72 °C, 10 min
4 °C

Once the PCR is finished, continue with the developing of the strip. If the samples are not to be tested immediately, they must be stored at 2-8 °C for no more than 24 hrs. For longer storage, freeze at -20 °C.

→ See the PCR preparation diagram in the Annex of the document.

2.- Strip developing

On request, Operon will provide the user with the working procedure for the specific instrument to be used for processing the strips (PST-60HL, Profiblot, Tendigo, Autoblot 3000H or Dynablot Heat). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>

Each PCR will be analyzed with 3 strips: strip A, strip B and strip C. Reagents and procedures are the same. For the denaturing process, **12.5 µl of denaturant + 12.5 µl** of PCR will be used per strip.

The product has been tested using different work platforms. For your information, the commercial names of the work platforms used are specified, which does not imply an exclusive use of such trademarks but the work protocol. The product works correctly as long as the platform used ensures the temperature program raised for the product.

→ See the Strip developing diagram in the Annex of the document.

3.- Strip interpretation

Strip interpretation can be done visually, using the evaluation chart included in the kit (see working procedures for each instrument) or automatically, with BLOTRix R2 or BLOTRix S1 scanners (BioSciTec).

On request, OPERON, S.A. will provide the user with the instructions for the automatic strip interpretation with the scanners (DO-09053002 "Instructions of use BLOTRix R2_S1"). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>

Important: results of strips interpretation with the scanners must always be visually verified; very faint bands could be not detected by the instruments.

RESULTS

S. PneumoStrip kit allows for the detection and genotyping of 76 of the most common serotypes of *S. pneumoniae*:


- Strip A: 1, 3, 4, 5, 6A, 6B, 6C, 6D, 7F/7A, 9A/9V, 14, 18A, 18B/18C, 18F, 19A, 19F, 23 F.
- Strip B: 2, 8, 9N/9L, 10A, 10B, 10F/C, 11A/D, 11B, 11C, 11F, 12F, 15B/15C, 17F, 20, 22F/22A, 33F/33A.
- Strip C: 7B, 7C/40, 15A, 15F, 16F, 19B/19C, 21, 25A/25F, 38, 24A, 24B/24F, 31, 32A/32F, 33B/33D, 33C, 35A, 35C, 35F, 47F, 41A, 41F.

A DNA sample may render the following results:


- No band: the bacteria is not identified as *S. pneumoniae*. It is important in these cases to verify that the PCR has worked properly (positive for control sample).
- Only the bands corresponding to the control genes: the bacteria is *S. pneumoniae* but the serotype does not correspond to any of those detected with the kit. In these cases, the bands corresponding to the genes *lyt* and *ply* will be detected, and the band associated with the CpSA gene may or may not be present (depending on whether or not it is an encapsulated or non-encapsulated bacterium). Attention: due to the presence of mutations and deletions that avoid the amplification of the CpSA gene, serotypes 25A, 25F and 38, independent of if they are capsulated or not, will not show the band associated to this gene (Leung M. et al. 2012).
- The control band and one/two bands associated to a certain type of *S. pneumoniae*: sample positive for that particular serotype. CpSA may or may not be present, depending on whether the sample is encapsulated or non-encapsulated (with the exception of serotypes 25A, 25F and 38, in which CpSA band will be always absent).
- The control band and several bands associated to different serotypes: positive sample with multiple serotypes.

Each serotype is identified with a band (aside from the control) except for:


- Serotypes from group 6: the identification of serotypes 6A, 6B, 6C and 6D is determined from a combination of three bands (6AC, 6BD and 6CD). If the sample corresponds to serotype 6A, the 6AC band will appear. If the sample corresponds to serotype 6B, the 6BD band will appear. If the sample corresponds to serotype 6C, the 6AC and 6CD bands will appear, and if the sample corresponds to serotype 6D, the 6BD and 6CD bands will appear.



- Serotypes from group 12: the identification of serotypes of group 12 is determined from a combination of three probes that allow to distinguish between 12A/46, 12B/44 and 12F.


- Serotypes from group 18: the identification of serotypes 18A, 18B/C and 18F is determined from a combination of two bands (the serotype 18 band common to all subtypes, and the band corresponding to each specific subtype). The presence of the first band (serotype 18) without the band associated to a subtype is indicative of the presence of a mutation in the probe binding zone. In these cases, the presence of serotype 18 is confirmed, but the subtype to which it belongs cannot be identified.

- Serotypes 15A and 15F; 24A, 24B and 24F; 35F and 47 F, 41A and 41F: the identification of these serotypes is determined from the combination of two bands, one band common to both serotypes and one band specific to one of them.

More serotyping examples using *S. PneumoStrip*

QUALITY CONTROL

The development control line should always appear (above the upper black line in strips A and B). Its absence will be indicative of issues encountered during the product hybridization/development stage.

The strip corresponding to the amplification negative control will show only the developing control band.

The strip corresponding to the positive amplification control will show the bands associated with the serotype present in the sample.

FEATURES / EXTERNAL AND INTERNAL EVALUATIONS

The use of this test is restricted to professional users that are familiar with the methods used in Molecular Biology.

The test has been evaluated with several DNA samples from different serotypes of S. pneumoniae.

- **Evaluation 1**

Conducted by Dr. José María Marimón in the Department of Microbiology of the Hospital Universitario de Donostia [University Hospital of Donostia] (San Sebastián).

129 samples of DNA were analysed from:

- 84 reference strains: DNA extraction from plate culture using the automated Nuclisens EasyMag extractor (BioMerieux).
- 45 samples from the routine operations of the microbiology laboratory of the Hospital Universitario de Donostia. DNA extraction from plate culture after boiling.

As a reference test, fragment analysis and the Quellung technique were used, according to the following procedure:

- a) All the samples were analysed with the S. PneumoStrip kit and with the fragment analysis serotyping method used by Dr. Marimón's group in routine diagnosis at the Hospital Universitario de Donostia. This method is based on the analysis of amplicons obtained from multiplexed capillary electrophoresis PCR, as described in Muñoz-Almagro C. *et al.*, 2014.
- b) The serotypes obtained with the multiplex PCR methods and S. PneumoStrip strips are verified using the Quellung technique.

The overall concordance obtained between the three serotyping methods was 99,2 % (128/129 samples). The serotyping concordance with the reference strains was 100 % and the serotyping concordance with the routine samples was 97,8 %, detecting 93.3% of serotypes present in this group of samples (reflecting the prevalence of the area of origin).

- **Evaluation 2**

Conducted at Operon S.A. (Zaragoza, Spain).

A total of 113 DNA samples were evaluated, previously analysed by fragment analysis (multiplex PCR + fragment analysis in a sequencer) and Quellung. For the study, both reference strains (CDC and SSI) and diagnostic strains (obtained from routine operations at the Hospital Universitario de Donostia) were used.

The reference strains were serotyped using the Quellung technique and the multiplex PCR system used routinely in the Department of Microbiology of the hospital.

The diagnostic strains were serotyped using the Quellung technique (groups 11, 10, 18) and/or the hospital's PCR multiplex system.

The percentage of concordance between methods was 100%.

POTENTIAL ISSUES

1. No bands appear on the strip, including the control band.

- The conjugate and/or developer were not balanced at room temperature.
- The conjugate and/or developer were not added, or were added in too small quantity.

2. Only the control band appears

- PCR amplification has failed (check with an agarose gel).
- The PCR was not added, or was added in too small quantity in the hybridization step.
- Incorrect hybridization or wash 1 buffers temperature (higher than instructed).
- Incorrect incubator temperature (higher than instructed)

3. The strip has a strong background colour following development

- The washing steps were not performed effectively (wrong timing, buffer quantity too small, cold buffers).

4. Non-homogenous strip staining

- Inadequate shaking (please review the plate thermo-shaker programme).
- The strips were not entirely submerged during the incubations.

5. Unexpected results

- Inadequate incubation and/or buffer/reagent temperatures.
- PCR contamination (verify using a negative control).
- Contamination of adjacent channels owing to the transfer of liquid from one channel to another when adding the hybridization buffer.

The working temperature range during the hybridisation process was evaluated with two thermoshakers (PST 60HL and Profiblot).

In the case of the Profiblot, the optimal working temperature range is between 42-43 °C.

The optimal working temperature range in the PST 60 HL is between 42-44 °C.

The probe most sensitive to temperature changes is the 6AC probe, for which reduction in working temperature could affect its specificity. The use of temperatures above those indicated could affect the sensitivity of the kit.

SENSITIVITY

We define **the limit of detection** of the kit as the minimum amount of DNA/RNA that can be detected with the test.

We define **the test limit** as the minimum amount of DNA/RNA that we can always genotype or always detect with the test.

Both limits have been determined for each one of the serotypes detected with the S. PneumoStrip kit. For this, 1/10 serial dilutions of DNA were analysed from each one of them, with three different product batches.

The sensitivity results obtained are shown in the following table:

Serotype	Limit of detection (ng/μl)	Serotype	Limit of detection (ng/μl)	S Serotype	Limit of detection (ng/μl)
1	0,0001	2	0,001	7B	0,0001
3	0,000001	8	0,0001	7C/40	0,0001
4	0,0001	9N/9L	0,00001	15A	0,00001
5	0,0001	10A	0,0001	15F	0,001
6A	0,0001	10B	0,0001	16F	0,0001
6B	0,0001	10F/C	0,0001	19B/19C	0,0001
6C	0,0001	11A/D	0,00001	21	0,0001
6D	0,0001	11B	0,0001	25A/25F	0,0001

Serotype	Limit of detection (ng/μl)	Serotype	Limit of detection (ng/μl)	S Serotype	Limit of detection (ng/μl)
7F/7A	0,0001	11C	0,0001	38	0,0001
9V/9A	0,0001	11F	0,00001	24A	0,00001
14	0,00001	12A/46	0,0001	24B/24F	0,0001
18A	0,0001	12B/12F/44	0,00001	31	0,0001
18B/18C	0,0001	12F	0,00001	32A/32F	0,0001
18F	0,00001	15B/15C	0,0001	33B/33D	0,0001
19A	0,0001	17F	0,00001	33C	0,0001
19F	0,00001	20	0,0001	35A	0,0001
23A	0,00001	22F/22A	0,00001	35C	0,0001
23B	0,0001	33F/33A	0,00001	35F	0,00001
23F	0,00001	37	0,00001	47F	0,00001
				41A	0,00001
				41F	0,00001

SPECIFICITY

The absence of reaction was checked with the following bacteria: *S. pseudoneumoniae*, *S. oralis*, *S. sanguis*, *S. pyogenes*, *S. agalactiae*, *S. anginosus*, *S. mitis*, *S. salivarius*, *E. coli* ST131, *E. coli* O6, *Klebsiella pneumoniae*, *Moraxella catharralis*, *Haemophilus influenzae*, *Legionella pneumophila*, *Staphylococcus aureus* (*mec A* + and *mec A* -).

HOOK EFFECT

The performance of *S. PneumoStrip* kit in the presence of increasing amounts of DNA (up to 100 ng) was evaluated. No PCR inhibition effect, or reduction in the intensity of the bands obtained, or problems arising from cross-reaction were observed. Therefore, there is no hook effect.

INTRA-ASSAY PRECISION

The intra-assay precision of *S. PneumoStrip* has been studied examining, first, 5 sensitivity curves (1/4 dilutions of a control DNA consisting of a mixture of DNA from all the serotypes detected by the kit) and on the other side, five replicates of 8 DNA samples positive for several serotypes.

Overall, the kit has shown a good intra-assay precision, obtaining the same serotyping for the samples and a similar sensitivity.

INTER-DAY PRECISION

The inter-day precision of *S. PneumoStrip* has been studied examining, first, a sensitivity curve (1/4 dilutions of a control DNA consisting of a mixture of DNA from all the serotypes detected by the kit) and on the other side, 8 real DNA samples. Each sensitivity curve was constructed by the same operator, using the same product lot and over five consecutive days.

Overall, the kit has shown a good inter-day precision, obtaining the same serotyping for the samples and a similar sensitivity.

INTER-LAB PRECISION

The inter-lab precision has been studied examining a sensitivity curve (1/4 dilutions of a control DNA consisting of a mixture of DNA from all the serotypes detected by the kit) and 8 real DNA samples. Each curve was constructed by a different operator, on the same day and using the same product lot. To vary the most of the conditions of use of the kit, each operator used a thermal cycler and different thermoshaker.

Overall, the kit has shown a good inter-lab precision, obtaining the same serotyping for the samples and a similar sensitivity.

INTER-LOT PRECISION

The inter-lot precision of S. PneumoStrip kit has been studied by the analysis using three different product lots of a sensitivity curve (1/4 dilutions of a control DNA consisting of a mixture of DNA from all the serotypes detected by the kit) and 25 DNA samples positive for some of the serotypes detected by the kit.

The kit has shown a good inter-lab precision. The results obtained for the 25 DNA samples were identical in all cases, and a similar sensitivity.

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1. Bentley S.D. et al. "Genetic analysis of the capsular biosynthetic locus from all 90 pneumococcal serotypes". PLOS Genetics (2006). Vol 2 (3): 262-269.
2. Hirst R.A. et al. "The role of pneumolysin in pneumococcal pneumonia and meningitis". Clin. Exp. Immunol. (2004). Vol 138: 195-201.
3. Huo Z, Spencer O, Miles J, Johnson J, Holliman R, Sheldon J, et al. Antibody response to pneumolysin and to pneumococcal capsular polysaccharide in healthy individuals and *Streptococcus pneumoniae* infected patients. Vaccine. 2004; 22:1157-61.
4. Leung M, Bryson K, Freystatter K, Pichon B, Edwards G, Charalambous B, Gillespie S. "Sequotyping: Serotyping *Streptococcus pneumoniae* by a single sequencing strategy". J.Clin.Microbiol (2012). Vol 50, N° 7: 2419-2427.
5. McAvin J. et al. "Sensitive and specific method for rapid identification of *Streptococcus pneumoniae* using real time PCR". J.Clin.Microbiol (2001). Vol 39, N° 10: 3446-3451.
6. Muñoz-Almagro C. et al. "Direct identification of *Streptococcus pneumoniae* capsular types in pleural fluids by using multiplex PCR combined with automated fluorescence-based capillary electrophoresis". Journal of Clinical Microbiology (2014), vol 52, n° 7, pags. 2736-2737.
7. Wasserman RL, Sorensen R. Evaluating children with respiratory tract infections: the role of immunization with bacterial polysaccharide vaccine. Pediatr Infect Dis J. 1999; 18:157-163.

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OPERON, S.A. - Camino del Plano, 19 - E-50410 Cuarte de Huerva - (Zaragoza) – SPAIN