IVD

High + Low PapillomaStrip

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EN

HIGH+LOW PAPILLOMASTRIP

Test for the detection of 37 anogenital HPVs

INTENDED USE

PapillomaStrip test is a test based on the reverse blot technique that allows qualitative detection in DNA samples from cervico-uterine smears or biopsies, of 37 human *Papillomavirus* subtypes:

- Medium-high risk: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 69, 73 and 82 (MM4 e IS39).
- Low risk: 6, 11, 40, 42, 43, 44, 54, 61, 62, 67, 70, 71, 72, 74, 81, 83, 84 y 91.

Take also under consideration other parameters such us symptoms or clinical history when giving the final diagnosis.

RATIONALE

High + Low PapillomaStrip kit is based on the principle of reverse hybridization and allows for the detection and identification of 37 genital HPVs in DNA samples from cervico-uterine smears or biopsies.

Human *Papillomavirus* (HPV) is responsible for a variety of skin and mucosa disorders. The clinical manifestations may vary according to the anatomic location of the lesion and the type of HPV involved. To date, over 100 different HPV genotypes have been described, of which more than 40 infect the anogenital tract. Of these, approximately 1/3 is associated with cervical cancer and anal neoplasm. The anogenital HPVs are generally divided into two categories: those with a low oncogenic potential risk (Low-Risk Group) and those with a medium-high oncogenic potential risk (Medium-High-Risk Group). The high-risk HPVs are generally associated with high-grade precancerous lesions and invasive cancer, while low-risk HPVs are frequently found in asymptomatic or benign conditions such as genital warts (1, 2).

The most common types of HPV are:

- Medium-high risk: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 69, 73 and 82.
- Low risk: 6, 11, 40, 42, 43, 44.

The HPV virus contains a single molecule of circular double stranded DNA of approximately 8 Kb. Three regions can be distinguished: a region regulating the transcription and replication processes (URR), a gene region of "early" expression encoding proteins involved in transcription, transformation and replication (E1-E2-E5-E6-E7) and a gene region of "late" expression encoding structural proteins (L1-L2-E4).

Currently, the diagnosis of HPV infection is based on viral DNA detection. The presence of conserved regions in different portions of viral DNA has allowed for the development of consensus primers, such as the systems MY09-MY11 (L1 region), PGMY09-PGMY11 (L1 region), Gp5+-Gp6+ (L1 region), SPF10 (L1 region) and LCR-E7 (E7 region), that enable the detection of a wide spectrum of different HPV genotypes. However, the differences in the oncogenic potential of HPVs make it necessary to identify the type of HPV present in the sample, for which different techniques are used: restriction analysis, sequencing, ELISA, reverse hybridization or microarrays (3, 4, 5, 6, 7, 8).

The procedure of High + Low PapillomaStrip kit consists in three steps: a) DNA extraction, b) amplification by PCR and c) hybridization/developing.

a) DNA extraction

Reagents not included in the kit.

Any standard DNA purification protocol can be used.

b) Amplification by PCR

High + Low PapillomaStrip kit includes all reagents required for the amplification of regions E6-E7 from the 37 HPVs detected by this kit. Since the expression of the E6-E7 genes is directly associated with carcinogenesis, this region has been selected as the detection of these targets may provide significant biological and pathological information. Furthermore, some studies show the existence of a gene deletion process in region L1. This deletion occurs when DNA from high-risk HPV integrates into the epithelial regions, and it has been described that it occurs in approximately 30% of the cases of cervical cancer samples positive for E6 region (9, 10).


Specific primers are used for each of the 37 HPVs in order to avoid problems that may occur with consensus primers for preferred amplification of any particular subtype.

During PCR, the amplification of a control gen (GADPH or β -globine) also occurs, acting as an indicator of the absence of PCR inhibitors in the amplified DNA sample.

c) Hybridization and developing

At this stage there is a specific binding of DNA fragments amplified during the PCR reaction to a series of probes deposited on nylon membranes.

Probes for each HPV detected are covalently linked to the membrane, along with a probe for the amplification control (GADPH or β -globine), a probe for the developing control and two or three lines, one black and one or two red, to control the positioning of the strip and as an aid in the interpretation of results (see attached diagram).

The detection of the fragments hybridized to different probes is performed using a conjugate (streptavidin-peroxidase) that binds to a biotin label that is added to the DNA fragments amplified during the PCR. Following the addition of a substrate for the peroxidase (TMB) a blue precipitate generates where hybridization has occurred.

As final result of the test, a band pattern is obtained which is interpreted with the aid of a control strip.

MATERIALS PROVIDED WITH THE KIT

HIGH + LOW PAPILOMASTRIP			16-test Kit	48-test Kit
	Membranes	STRIPS	16 (8 + 8)	48 (24 + 24)
	8-channel trays	PL	2	-
PCR Reagents REAG PCR	PCR Mixture	REAG PCR MIX	1 ml	1.2 x 2 ml
	Primers	REAG PRIMER	0.07 ml 0.07 ml	0.15 ml 0.15 ml
	Taq	REAG TAQ HS	0.05 ml	0.12 ml
	Denaturant	SOLN DN	1 ml	1 ml
	Hybridization Buffer	BUF HYB	60 ml	125 ml
	Wash Buffer 1	BUF WASH 1	100 ml	250 ml
	Conjugate	CONJ HRP	60 ml	125 ml
	Wash Buffer 2	BUF WASH 2	130 ml	350 ml
	Substrate	SUBS TMB	30 ml	65 ml
	Instructions of use		1	1
	Interpretation Chart		1 + 1	1 + 1

MATERIALS NOT INCLUDED IN THE KIT

The following additional material is required when using the kit:

1. Microtubes for PCR
2. Micropipettes and micropipette tips (sterile or UV-irradiated and ideally with a filter)
3. Tweezers and a pencil (optional)
4. Chronometer
5. Thermometer

REQUIRED EQUIPMENTS FOR KIT DEVELOPMENT

1. Thermoshakers for plates
The kit has been validated with the following thermoshakers: PST-60HL (Biosan S.L.), Profiblot T30 y T48 (Tecan), Autoblot 3000H (medTEC), TENDIGO (Fujirebio), Dynablot Heat (Dy nex).
2. Thermocycler
The following thermocyclers have been successfully used with the kit: PTC-100 (MJ Research), MJ Mini Gradient Thermal Cycler (BioRad), Mastercycler Personal (Eppendorf), S1000 Thermal Cycler (BioRad).
3. BLOTri x S1 or BLOTri x R2 scanner (BioSciTec, optional).

PRECAUTIONS

1. Only use the reagents *in vitro*.
2. Strictly follow the instructions provided. Modifying any of the prescribed **steps or temperatures** may severely affect the test results.
3. All the materials and reagents used must be free of DNAses. The use of filtered pipette tips is recommended for PCR preparation to prevent aerosol contamination. Ensure that all the usual precautions are taken in the preparation of the amplifications. Use autoclaved material for the hybridization/development process.
4. Store the kit components as indicated in the instructions.
5. Do not exchange components from kits with different lot numbers.

6. Do not use kit components after the expiration dates.
7. Strips are for single-use.
8. Supplied trays are for single-use.
9. Each channel of the supplied trays is for just one strip.
10. If the package is broken, the product can still be used providing none of the components have been damaged.
11. The used product should be discarded in compliance with current legislation.
12. Patient samples must always be treated as potentially infectious. Environmental and safety standards must be adhered to.
13. Once it is used, it is advisable to manipulate the strip taking into account the same considerations as for the sample. The strips have to be manipulated with the proper precautions and should be managed as biohazardous material.
14. The denaturation solution contains < 2% NaOH and is an irritant to both eyes and skin (H314 and P280, P305, P351, P338, P310).
15. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking and component lots.

STORAGE

Store all reagents at 2-8 °C. Precipitate formation may occur in some solutions (hybridization buffer, wash buffer 1 and 2, conjugate) due to storage at 2-8 °C. These precipitates are redissolved when brought to room temperature (conjugate) or heated to 42 °C hybridization buffer and wash buffer 1 and 2).

Expiry date for all reagents is printed on the respective label.

SAMPLES

The test has been designed and validated for its use with DNA derived from cervico-uterine smears. It has also shown to perform correctly with samples from penis swabs. The use of samples from another source must be validated by the user

For sampling use dry and sterile cotton swabs or brushes. Make sure to take sufficient sample amount, but without causing bleeding from the lesion. Put the swab into its tube, do not use any preserving medium. Store at 2-8 °C if it is to be processed within one week or at -20 °C for a later processing.

The quality and concentration of DNA extracted from cervico-uterine specimens is crucial for the proper functioning of the test. For this reason, the kit includes a control that detects the presence of PCR inhibitors in DNA (amplification of GAPDH or β -globine genes). Additionally, we recommend measuring the concentration of DNA samples to be used, since the use of DNA samples with a concentration below 10 ng/ μ l, although showing control gene amplification, could reduce the sensitivity of the test.

Some of the DNA extraction procedures that have been evaluated with satisfactory results with PapillomaStrip kit are:

- a) DNA extraction kit from OPERON S.A.: Ref 3.149.020.64.000, DNA extraction kit.
- b) DNA extraction with silica resins (QIAamp DNA Mini Kit; Qiagen).
- c) Extraction protocol for microarray Clart HPV2 from Genomica S.L.
- d) DNA extraction with magnetic particles.

Store the DNA samples at 2-8 °C if they are going to be used shortly after, or at -20 °C for longer periods of storage.

HIGH + LOW PAPILOMASTRIP PROCEDURE

1.-Polymerase chain reaction

Important: before opening the vials with the PCR reagents centrifuge them briefly. This will ensure that all the contents will be at the bottom of the tube.

PCR preparation

For each sample two amplification reactions will be run: one of them with medium-high risk primers (vial with blue cap) and another one with low risk primers (vial with white cap). PCR Premix and Taq reagents are common to both reactions.

- Place the tubes required according to the number of DNAs to be amplified (two per sample).
- Add to each PCR tube: 38 μ l of PCR premix + 5 μ l of primers (medium-high or low risk) + 2 μ l of Taq + 5 μ l of DNA. Mix. Whenever possible, keep all reagents at 2-8 °C during preparation. If several DNAs are going to be amplified, it is recommended to prepare a common mixture containing all reagents, so 45 μ l premix + 5 μ l DNA are finally added. For example:

Nr of PCRs	PCR Premix	Primers	Taq
1	38 μ l	5 μ l	2 μ l
3	133 μ l	17.5 μ l	7 μ l
5	209 μ l	27.5 μ l	11 μ l
8	342 μ l	45 μ l	18 μ l

* Always prepare an excess of the mixtures containing all PCR reagents, to compensate volume losses during pipetting process.

Due to requirements of "good laboratory practice", the user should also include a negative control (e.g. water or TE buffer) to exclude contamination and, if required, a positive control (not included in the kit).

Amplification

Regard the importance of getting a suitable concentration of the DNA to get a proper result.

Place the tubes into the cycler device (if required, add 1 drop of Nujol oil on top of every PCR reaction) and amplify DNA by the following program:

94 °C , 5 min
40 cycles of: 94 °C, 1 min
58 °C, 1 min
72 °C, 1 min
72 °C, 5 min
4 °C

Once the PCR is finished, continue with the developing of the strip. If the samples are not to be tested immediately, they must be stored at 2-8 °C for no more than 24 hrs. For longer storage, freeze at -20 °C.

→ See the PCR preparation diagram in the Annex of the document.

2.- Strip developing

The analysis of the PCRs amplified with low-medium risk primers is conducted using the HIGH PapillomaStrip strip. That of the PCRs amplified with the low risk primers is done with the LOW PapillomaStrip strips. In both cases, the reagents and procedures to follow are common. For the denaturing process, 12.5 μ l of denaturant + 12.5 μ l of PCR will be used.

On request, Operon will provide the user with the working procedure for the specific instrument to be used for processing the strips (PST-60HL, Profiblot, Autoblots 3000H, TENDIGO or Dynablot Heat). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>

The product has been tested using different work platforms. For your information, the commercial names of the work platforms used are specified, which does not imply an exclusive use of such trademarks but the work protocol. The product works correctly as long as the platform used ensures the temperature program raised for the product.

→ See the Strip developing diagram in the Annex of the document.

3.- Strip interpretation

Strip interpretation can be done visually, using the evaluation chart included in the kit (see working procedures for each instrument) or automatically, with BLOTRix R2 or BLOTRix S1scanners (BioSciTec).

On request, OPERON, S.A. will provide the user with the instructions for the automatic strip interpretation with the scanners (DO-09053004 "Instructions of use BLOTRix R2_S1"). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>

Important: results of strips interpretation with the scanners must always be visually verified; very faint bands could be not detected by the instruments.

RESULTS

High + Low PapillomaStrip kit allows for the detection and genotyping of 18 low risk anogenital HPVs and 19 high-medium risk HPVs:


- High-medium risk: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 69, 73 y 82 (MM4 e IS39).
- Low risk: 6, 11, 40, 42, 43, 44, 54, 61, 62, 67, 70, 71, 72, 74, 81, 83, 84 y 91.

The bands obtained on each strip are interpreted with the aid of two charts: one of them for HPVs low risk strip and another one for HPVs high risk strip.

A DNA sample may render the following results:

- The control band only (GADPH for low risk, or β -globin for medium-high risk): none of the 37 HPVs detected by this kit is present or could be amplified.
- The control band and one band associated to a certain type of HPV: sample positive for that particular HPV.
- The control band and several bands associated to different HPVs: positive sample with multiple infection.
- One or several bands associated to different HPVs: positive sample with infection by one or several HPVs in which the GADPH or β -globine genes do not amplify. This may happen, for example, when the sample has a high number of copies of one or more viruses, as the PCR conditions have been designed to favor the amplification of HPV DNA versus that of the internal control.

Examples of genotyping using High + Low PapillomaStrip

QUALITY CONTROL

The development control line should always appear above the upper black line. Its absence will be indicative of issues encountered during the product hybridization/development stage.

The strip corresponding to the amplification negative control will show only the developing control band. For negative samples, the presence of the band associated to the internal control (GADPH or β -globin) is required to confirm that it is a true negative. The absence of this band in a negative sample is indicative of the absence of amplification during the PCR, which can be caused by lack of DNA (measure concentration), by degradation or by the presence of PCR inhibitors in the same (hemoglobin, traces of salts, etc.).

FEATURES / EXTERNAL AND INTERNAL EVALUATIONS

The use of this test is restricted to professional users that are familiar with the methods used in Molecular Biology.

The test has been evaluated with several DNA samples from cervico-uterine smears, obtained by different protocols, with satisfactory results. The test has been successfully used for the analysis of the presence of HPV in penis swabs.

Evaluation nº 1

The sensitivity and specificity of High + Low PapillomaStrip were compared with the Human Papillomavirus CLART test from Genomica (Zeltia Group). The study was performed at Operon and the Centro de Análisis Genéticos C.A.G.T. (Zaragoza). A total of 90 samples of cervico-vaginal smears were analyzed. The DNA was extracted from this cervico-vaginal smears according to the protocol of the Human Papillomavirus CLART kit.

The HPVs detected with the Human Papillomavirus CLART kit were:

- Low risk: 6, 11, 40, 42, 44, 54, 61, 62, 71, 72, 81, 83, 84.
- High risk: 16, 18, 26, 31, 33, 35, 39, 43, 45, 51, 52, 53, 56, 58, 59, 66, 68, 70, 73, 82, 85, 89.

The results are shown in the table below:

Number of analyzed samples: 90	
Negative for microarray: 32	Negative Low PapillomaStrip: 27 (84.4 %) Positive Low PapillomaStrip: 5 (15.6 %) Negative High PapillomaStrip: 27 (84.3 %) Positive High PapillomaStrip: 5 (15.6 %) Negative High+Low: 24 (75 %)
Positive for microarray – High/medium risk HPV: 34	Negative Low PapillomaStrip: 15 (44,1 %) Positive Low PapillomaStrip: 19 (55.9 %) Negative High PapillomaStrip: 1 (2.9 %) Positive High PapillomaStrip: 33 (97.1 %) Negative High+Low: 1 (2.9 %)
Positive for microarray – Low risk HPV: 7	Negative Low PapillomaStrip: 0 (0 %) Positive Low PapillomaStrip: 7 (100 %) Negative High PapillomaStrip: 4 (57.1 %) Positive High PapillomaStrip: 3 (42.9 %) Negative High + Low: 0 (0 %)
Positive for microarray – High/medium risk and/or low risk: 17	Negative Low PapillomaStrip: 0 (0 %) Positive Low PapillomaStrip: 17 (100 %) Negative High PapillomaStrip: 0 (0 %) Positive High PapillomaStrip: 17 (100 %) Negative High + Low: 0 (0 %)

	N° of positive samples			% of detection in sample		
	High Pap	CLART	Total	High Pap	CLART	Total
16	15	11	15	16.7	12.2	16.7
18	3	1	3	3.3	1,1	3.3
26	0	0	0	0	0	0
31	14	10	14	15.6	11.1	15.6
33	4	4	5	4.4	4.4	5.6
35	5	5	6	5.6	5.6	6.7
39	10	4	10	11.1	4.4	11.1
45	1	1	1	1.1	1.1	1.1
51	16	15	16	17.8	16.7	17.8
52	9	5	9	10	5.6	10
53	12	12	13	13.3	13.3	14.4
56	5	1	5	7.6	1.5	7.6
58	2	3	3	3	4.5	4.5
59	3	3	3	4,5	4,5	4,5
66	12	14	16	13.3	15.6	17.8
68	6	0	6	6.7	0	6.7
69	0	0	0	0	0	0
73	6	0	6	6.7	0	6.7
82	2	2	2	2.2	2.2	2.2

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Evaluation n° 2

The sensitivity and specificity of High + Low PapillomaStrip were compared with the InnoLipa HPV Genotyping test (Innogenetics). The study was performed at OPERON (High + Low PapillomaStrip) and Genyca Innova Laboratories (InnoLipa HPV Genotyping, Madrid).

The HPVs detected with the Innogenetics kit were:

- Low risk: 6, 11, 40, 43, 44, 54, 70, 71, 74.
- High risk: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68/73, 69, 82.

The results are shown in the table below:

Number of analyzed samples: 40	
Positive samples in InnoLipa – Low risk HPV: 15	Negative Low PapillomaStrip: 0 (0 %) Positive Low PapillomaStrip: 15 (100 %) Negative High PapillomaStrip: 11 (73,3 %) Positive High PapillomaStrip: 4 (26,7 %) Negative High + Low: 0 (0 %)
Positive samples in InnoLipa – High/Medium risk and Low risk HPVs: 25	Negative Low PapillomaStrip: 2 (8 %) Positive Low PapillomaStrip: 23 (92 %) Negative High PapillomaStrip: 1 (4 %) Positive High PapillomaStrip: 24 (96 %) Negative High + Low: 0 (0 %)
<p>Interpretation of results</p> <p>In 26.7% of samples in which the strip of Innogenetics only detected low-risk HPVs, with the kit High PapillomaStrip we detected any high-risk HPV.</p> <p>Among the samples with high and low risk HPVs, we did not detect two low-risk HPVs and one high risk HPV. In 100% of these samples some HPV is detected.</p>	

	N° of positive samples			% of detection in sample		
	High Pap	Inno	Total	High Pap	Inno	Total
16	4	3	6	10	7,5	15
18	3	1	3	7,5	2,5	7,5
26	2	1	2	5	2,5	5
31	3	2	3	7,5	5	7,5
33	3	3	3	7,5	7,5	7,5
35	3	3	3	7,5	7,5	7,5
39	2	0	2	5	0	5
45	3	0	3	7,5	7,5	7,5
51	9	9	9	22,5	22,5	22,5
52	4	3	4	10	7,5	10
53	4	5	5	10	12,5	12,5
56	4	2	4	10	5	10
58	4	1	4	10	2,5	10

	N° of positive samples			% of detection in sample		
	High Pap	Inno	Total	High Pap	Inno	Total
59	2	0	2	5	0	5
66	4	3	4	10	7,5	10
68	1	0	1	2,5	0	2,5
69	0	0	0	0	0	0
73	0	0	0	0	0	0
82	0	0	0	0	0	0

	N° of positive samples			% of detection in sample		
	Low Pap	Inno	Total	Low Pap	Inno	Total
6	16	16	16	40	40	40
11	5	6	6	12,5	15	15
40	0	0	0	0	0	0
42	10	NA	10	25	NA	25
43	4	3	5	10	7,5	12,5
44	7	8	8	17,5	20	20
54	9	7	9	22,5	17,5	22,5
61	2	NA	2	5	NA	5
62	1	NA	1	2,5	NA	2,5
67	2	NA	2	5	NA	5
70	0	0	0	0	0	0
71	0	0	0	0	0	0
72	0	NA	0	0	NA	0
74	5	4	6	12,5	10	15
81	0	NA	0	0	NA	0
83	0	NA	0	0	NA	0
84	2	NA	2	5	NA	5
91	4	NA	4	10	NA	10

POTENTIAL ISSUES

1. No bands appear on the strip, including the control band.

- The conjugate and/or developer were not balanced at room temperature.
- The conjugate and/or developer were not added, or were added in too small quantity.

2. Only the control band appears

- PCR amplification has failed (check with an agarose gel).
- The PCR was not added, or was added in too small quantity in the hybridization step.
- Incorrect hybridization or wash 1 buffers temperature (higher than instructed).
- Incorrect incubator temperature (higher than instructed).

3. The strip has a strong background colour following development

- The washing steps were not performed effectively (wrong timing, buffer quantity too small, cold buffers).

4. Non-homogenous strip staining

- Inadequate shaking (please review the plate thermo-shaker programme).
- The strips were not entirely submerged during the incubations.

5. Unexpected results

- Inadequate incubation and/or buffer/reagent temperatures.
- PCR contamination (verify using a negative control).
- Contamination of adjacent channels owing to the transfer of liquid from one channel to another when adding the hybridization buffer.

The following table displays in gray the working temperature ranges (hybridization/developing; test performed using the PST-60HL thermoshaker) with which adequate results are obtained for each band (good hybridization intensity, no cross-reaction):

	39 °C	40 °C	41 °C	42 °C	43 °C	44 °C	45 °C		39 °C	40 °C	41 °C	42 °C	43 °C	44 °C	45 °C
SBG								S52							
S16								S53							
S18								S56							
S26								S58							
S31								S59							
S33								S66							
S35								S68							
S39								S69							
S45								S73							
S51								S82							

	39 °C	40 °C	41 °C	42 °C	43 °C	44 °C	45 °C		39 °C	40 °C	41 °C	42 °C	43 °C	44 °C	45 °C
GADPH								S67							
S6								S70							
S11								S71							
S40								S72							
S42								S74							
S43								S81							
S44								S83							
S54								S84							
S61								S91							
S62															

S11 y S74: detection sensibility is reduced at temperatures higher than 43 °C.
 S91, S43, S44: band intensity is reduced at temperatures higher than 43 °C, without affecting detection sensibility.

SENSITIVITY

The detection sensitivity of High + Low PapillomaStrip for each of the HPVs detected by the kit was determined.

For this purpose, plasmids of regions E6-E7 from the 37 HPVs were prepared in the pGemT vector. The concentration of each plasmid was quantified by spectrophotometry at 260 nm and the High + Low

PapillomaStrip detection sensitivity for each HPV was determined by the analysis of 1/10 dilutions of the different plasmids in a negative control DNA.

Sensitivity results obtained are shown in the following table:

HPV subtype	Sensitivity	
	Concentration (pg/ul)	pg per PCR
HPV16	0,0002	0,001
HPV18	0,0002	0,001
HPV26	0,00002	0,0001
HPV31	0,00002	0,0001
HPV33	0,00002	0,0001
HPV35	0,00002	0,0001
HPV39	0,0002	0,001
HPV45	0,00002	0,0001
HPV51	0,0002	0,001
HPV52	0,00002	0,0001

HPV subtype	Sensitivity	
	Concentration (pg/ul)	pg per PCR
HPV53	0,00002	0,0001
HPV56	0,00002	0,0001
HPV58	0,0002	0,001
HPV59	0,0002	0,001
HPV66	0,0002	0,001
HPV68	0,0002	0,001
HPV69	0,0002	0,001
HPV73	0,00002	0,0001
HPV82	0,002	0,01

HPV subtype	Sensitivity	
	Concentration (pg/ul)	pg per PCR
HPV6	0,0002	0,001
HPV11	0,0002	0,001
HPV40	0,00002	0,0001
HPV42	0,00002	0,0001
HPV43	0,0002	0,001
HPV44	0,00002	0,0001
HPV54	0,00002	0,0001
HPV61	0,00002	0,0001
HPV62	0,0002	0,001

HPV subtype	Sensitivity	
	Concentration (pg/ul)	pg per PCR
HPV67	0,00002	0,0001
HPV70	0,002	0,01
HPV71	0,00002	0,0001
HPV72	0,00002	0,0001
HPV74	0,00002	0,0001
HPV81	0,00002	0,0001
HPV83	0,00002	0,0001
HPV84	0,0002	0,001
HPV91	0,00002	0,0001

SPECIFICITY

No cross-reaction has been detected with the following pathogens agents:

<i>Escherichia coli</i>	<i>Streptococcus agalactiae</i>
<i>Neisseria gonorrhoeae</i> tipo I	<i>Staphylococcus aureus</i> mecA (-)
<i>Neisseria gonorrhoeae</i> tipo II/III	<i>Staphylococcus aureus</i> mecA (+)
<i>Herpes simplex</i> 1	<i>Cándida albicans</i>
<i>Herpes simplex</i> 2	<i>Ureaplasma urealyticum</i>
<i>Haemophilus ducreyi</i>	<i>Mycoplasma genitalium</i>
<i>Trichomonas vaginalis</i>	

HOOK EFFECT

The performance of the High + Low PapillomaStrip kit in the presence of increasing amounts of plasmidic DNA (up to 100 ng) was evaluated. No PCR inhibition effect, or reduction in the intensity of the bands obtained, or problems arising from cross-reaction were observed. Therefore, there is no Hook effect.

INTRA-ASSAY PRECISION

The intra-assay precision of High + Low PapillomaStrip has been studied examining, first, 5 sensitivity curves (1/10 dilutions of a control DNA consisting of a mixture of plasmids of the 18 low risk or the 19 high-medium risk HPVs detected by the kit) and on the other side, five replicates of five DNA samples obtained from cervico-vaginal smears.

Overall, High + Low PapillomaStrip kit has shown a good intra-assay precision. Both the test performed with dilution curves (where all detected HPVs are amplified simultaneously) and the test using actual samples, show a very good reproducibility for samples of medium, normal or strong intensity and a slight variation in samples near the limit of detection of the different HPVs.

INTER-DAY PRECISION

The inter-day precision of High + Low PapillomaStrip has been studied examining 5 sensitivity curves (1/10 dilutions of a control DNA consisting of a mixture of plasmids of the 18 low risk or the 19 high-medium risk HPVs detected by the kit). Each sensitivity curve was constructed by the same operator, using the same product lot and over five consecutive days.

Overall, the High + Low PapillomaStrip kit shows a good reproducibility varying the day. As in the intra-assay precision, a very good reproducibility for samples of medium, normal or strong intensity and a slight variability for the samples near the limit of detection of the different HPVs can be observed.

INTER-LAB PRECISION

The inter-lab precision has been studied examining 3 sensitivity curves (1/10 dilutions of a control DNA consisting of a mixture of plasmids of the 18 low risk or the 19 high-medium risk HPVs detected by the kit). Each curve was constructed by a different operator, on the same day and using the same product lot.

Overall, the High + Low PapillomaStrip kit shows a good reproducibility varying the operator and the conditions. The largest differences are found, as for the rest of the tests performed with these samples (mixtures of all detected HPVs), in the dilutions near the limit of sensitivity of each HPV.

INTER-LOT PRECISION

The inter-lot precision of High + Low PapillomaStrip kit has been studied by the analysis using three different product lots of a sensitivity curve (1/10 dilutions of a control DNA consisting of a mixture of plasmids of the 18 low risk or the 19 high-medium risk HPVs detected by the kit) and 37 plasmid samples from each of the HPVs detected by this kit.

Overall, the High + Low PapillomaStrip kit has shown a good inter-lot precision. The results obtained for the 37 plasmid samples were identical in all cases; the results of the sensitivity curves, as in the rest of the precision tests, have demonstrated a very good reproducibility for samples of medium, normal or strong intensity and a light variability for the samples near the limit of detection of the different HPVs.

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1. Bosch, F.X. et al. "Prevalence of human papillomavirus in cervical cancer: a worldwide perspective. International biological study on cervical cancer (IBSCC) Study Group". *J. Natl. Cancer Inst.* (1995). 87: 796-802.
2. Bosch, F.X. et al. "Papillomavirus research update: highlights of the Barcelona HPV 2000 international papillomavirus conference". *J. Clin. Pathol.* (2001). 54: 163-175.
3. Gravitt P.E. et al. "Genotyping of 27 human papillomavirus types by using L1 consensus PCR products by a single-hybridization, reverse line blot detection method". *J. Clin. Microbiol.* (1998). Vol 36, N° 10: 3020-3027.
4. Coutlèe F. et al. "Enhanced detection and typing of human papillomavirus (HPV) DNA in anogenital samples with PGM1 primers and the linear array HPV genotyping test". *J. Clin. Microbiol.* (2006). Vol 44, N° 6: 1998-2006.
5. Jacobs M.V. et al. "Group-specific differentiation between high- and low-risk human papillomavirus genotypes by general primer-mediated PCR and two cocktails of oligonucleotide probes". *J. Clin. Microbiol.* Vol 33, N° 4: 901-905.
6. Kleter B. et al. "Development and clinical evaluation of a highly sensitive PCR-reverse hybridization line probe assay for detection and identification of anogenital human papillomavirus". *J. Clin. Microbiol.* (1999). Vol 37, N° 8: 2508-2517.
7. Gheit T. "Development of a sensitive and specific assay combining multiplex PCR and DNA microarray primer extension to detect high-risk mucosal human papillomavirus types". *J. Clin. Microbiol.* (2006). Vol 44, N° 6: 2025-2031.
8. Sandri M.T. "Comparison of the digene HC2 assay and the Roche AMPLICOR human papillomavirus (HPV) test for detection of high-risk HPV genotypes in cervical samples". *J. Clin. Microbiol.* (2006). Vol 44, N° 6: 2141-2146.
9. Han, J. et al. "Simultaneous amplification and identification of 25 human papillomavirus types with Templex Technology". (2006). *J. Clin. Microbiol.* Vol 44, N°11: 4157-4162.
10. Sotlar K. et al. "Detection and typing of human papillomavirus by E6 nested multiplex PCR". (2004). *J. Clin. Microbiol.* Vol 42, N° 7: 3176-3184.

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OPERON, S.A. - Camino del Plano, 19 - E-50410 Cuarte de Huerva - (Zaragoza) – SPAIN