IVD

ThromboStrip

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EN

ThromboStrip

Test for the detection of point mutations associated with venous thrombosis risk

INTENDED USE

ThromboStrip test is a test designed for the detection of three point mutations/genetic variants associated with the risk of suffering deep venous thrombosis: G1691A mutation in Factor V; G20210A mutation in prothrombin; and C677T polymorphism in methylenetetrahydrofolate reductase.

This test is only a tool to help the professionals in the diagnosis of the indicated pathologies. This kit cannot be the only diagnosis tool.


BASIS

ThromboStrip kit is based on reverse hybridization principle and allows for the identification of three genetic variants (G1691A mutation in Factor V gene and G20210A mutation in prothrombin gene, and C677T polymorphism of MTHFR gene).

The procedure consists of three steps: a) DNA extraction (reagents not included in the kit), b) PCR amplification, and c) hybridization/detection.

Fragments of coagulation factor V (FV) gene, coagulation factor II (prothrombin, PTB) gene, and methylenetetrahydrofolate reductase (MTHFR) gene, containing G1691A, G20210A, and C677T mutations/polymorphism respectively, are simultaneously amplified by means of polymerase chain reaction. During this PCR, a label (biotin) is added to the amplified DNA fragments.

Test membranes bear covalently-linked DNA probes which specifically recognize each gene sequence amplified during the PCR. There are two probes for each gene: one for the normal sequence and one for the mutated sequence. A control probe of the developing process and a black and a red lines for the control of strip positioning and as a result interpretation aid are also included.

During hybridization process, the probes specifically bind to DNA fragments amplified during PCR.

Once PCR products are hybridized with membrane probes (and after several washes to eliminate any non-specific binding) the hybridization is detected ("developing"). With this aim the membrane is

incubated with a streptavidin-peroxidase conjugate, which will bind to every site with a biotin molecule; then, a peroxidase substrate (TMB) is added which will develop a blue precipitate in the site where hybridization has occurred.

As final result of the test, a band pattern is obtained which is interpreted with the aid of a control strip.

Each gene may present three situations (according to presence or not of the mutation, and, if present whether in homozygous or heterozygous status): N/N, N/M, M/M.

MATERIALS INCLUDED IN THE KIT

THROMBOSTRIP			16 tests kit
	Membranes	STRIPS	16
	8-channel plates	PL	2
PCR Reagents REAG PCR	PCR Mix	REAG PCR MIX	0.8 ml
	Primers	REAG PRIMER	0.11 ml
	Taq	REAG TAQ	0.030 ml
	Denaturant	SOLN DN	1 ml
	Hybridization buffer	BUF HYB	60 ml
	Wash buffer 1	BUF WASH 1	100 ml
	Conjugate	CONJ HRP	60 ml
	Wash buffer 2	BUF WASH 2	130 ml
	Substrate	SUBS TMB	30 ml
	Instructions of use		1
	Interpretation chart		1

MATERIALS NOT INCLUDED IN THE KIT

The following additional material is required when using the kit:

1. Microtubes for PCR.
2. Micropipettes and micropipette tips (sterile or UV-irradiated and ideally with a filter).
3. Tweezers and a pencil (optional).
4. Chronometer.
5. Thermometer.

REQUIRED EQUIPMENTS FOR KIT DEVELOPMENT

1. Thermoshakers for plates
The kit has been validated with the following thermoshakers: PST-60HL (Biosan S.L.), Profiblot T30 y T48 (Tecan), Autoblot 3000H (medTEC), TENDIGO (Fujirebio), Dynablot Heat (Dynex).
2. Thermocycler
The following thermocyclers have been successfully used with the kit: PTC-100 (MJ Research), MJ Mini Gradient Thermal Cycler (BioRad), Mastercycler Personal (Eppendorf) and S1000 Thermal Cycler (BioRad).
3. BLOTrix S1 or BLOTrix R2 scanner (BioSciTec, optional).

PRECAUTIONS

1. Only use the reagents *in vitro*.
2. Strictly follow the instructions provided. Modifying any of the prescribed **steps or temperatures** may severely affect the test results.
3. All the materials and reagents used must be free of DNAses. The use of filtered pipette tips is recommended for PCR preparation to prevent aerosol contamination. Ensure that all the usual precautions are taken in the preparation of the amplifications. Use autoclaved material for the hybridization/development process.
4. Store the kit components as indicated in the instructions.
5. Do not exchange components from kits with different lot numbers.
6. Do not use kit components after the expiration dates.
7. Strips are for single-use.
8. Supplied trays are for single-use.
9. Each channel of the supplied trays is for just one strip.
10. If the package is broken, the product can still be used providing none of the components have been damaged.
11. The used product should be discarded in compliance with current legislation.
12. Patient samples must always be treated as potentially infectious. Environmental and safety standards must be adhered to.
13. Once it is used, it is advisable to manipulate the strip taking into account the same considerations as for the sample. The strips have to be manipulated with the proper precautions and should be managed as biohazardous material.
14. The denaturation solution contains < 2% NaOH and is an irritant to both eyes and skin (H314 and P280, P305, P351, P338, P310).
15. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking and component lots.

STORAGE

Store at 2-8°C. Precipitate formation may occur in some solutions (hybridization buffer, wash buffer 1 and 2, conjugate) due to storage at 2-8 °C. These precipitates will redissolve during temper (conjugate) or warming to 42 °C (hybridization and wash buffers).

Expiry date for all reagents is printed on the respective label.

SAMPLES

The base material for carrying out the test is DNA of a known concentration. This DNA can be extracted from blood (with EDTA or citrate), saliva or tissue and should be of sufficient quality for PCR amplification (R 260/280 >1,5).

DNA samples, obtained by several DNA extraction methods, have been tested on ThromboStrip. The extraction procedures were the following:

- a) Salting out (DNA extraction by saline precipitation).
- b) DNA extraction with silica and silica columns.
- c) DNA extraction of a drop of blood in FTA paper with Chelex (8 mm diameter circle extracted with 25 µl of 10 % Chelex solutions; 5 µl of DNA per PCR).
- d) DNA extraction of a drop of blood in FTA paper with Whatmann solutions (2 mm diameter circle per PCR).

Store the DNA samples at 2-8 °C if they are going to be used shortly after, or at -20 °C for longer periods of storage.

THROMBOSTRIP PROCEDURE

1.-Polymerase chain reaction

PCR preparation

Important: before opening the vials with the PCR reagents centrifuge them briefly. This will ensure that all the contents will be at the bottom of the tube.

- Prepare the required PCR tubes based on the number of DNAs to be amplified.
- Add to each PCR tube: 39 μl of PCR premix + 5 μl of primers + 1 μl of Taq + 5 μl of DNA (50 ng/ μl). Mix well. If possible, keep the reagents and the mixtures at 2-8 °C during the preparation.

If various DNA samples are to be amplified, preparing a common mixture with all the reagents is recommended to finally add 45 μl of mixture + 5 μl of DNA (50 ng/ μl). For example:

Nr. of PCRs	PCR premix	Primers	Taq
1	39 μl	5 μl	1 μl
3	136.5 μl	17.5 μl	3.5 μl
5	214.5 μl	27.5 μl	5.5 μl
8	351 μl	45 μl	9 μl

* Always prepare an excess of the mixtures containing all PCR reagents, to compensate volume loss during pipetting process.

Due to requirements of "good laboratory practice", the user should also include a negative control (e.g. water or TE buffer) to exclude contamination and, if required, a positive control (not included in the kit).

Amplification

Regard the importance of getting a suitable concentration of the DNA to get a proper result.

Place the microtubes into the thermocycler (if the apparatus requires the use of Nujol oil, add 1 drop on top of every PCR reaction) and amplify DNA by following program:

94 °C, 3 min
30 cycles of: 94 °C, 1 min
54 °C, 1 min
72 °C, 1 min
72 °C, 5 min
4 °C

3. Verification of PCR

If desired, the result of the PCR can be verified in a 3 % agarose gel. For this, apply 5 μl of each PCR reaction in each gel well (PCR premix includes loading buffer and two dyes for electrophoresis control). Run the electrophoresis at 100 V for 1 hr and examine under UV light. Three bands of the following sizes should be detected: 245 bp, 234 bp and 206 bp (according to the size of the gel used, 245 and 234 bp bands may appear as a single band).

Once the PCR is finished, continue with the developing of the strip. If the samples are not to be tested immediately, they must be stored at 2-8 °C for no more than 24 hrs. For longer storage, freeze at -20 °C.

→ See the PCR preparation diagram in the Annex of the document.

2.- Strip development

On request, Operon will provide the user with the working procedure for the specific instrument to be used for processing the strips (PST-60HL, Profiblot, Autoblot 3000H, TENDIGO or Dynablot Heat). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>

In the denaturing step, **12.5 µl of denaturation buffer + 12.5 µl of PCR** will be dispensed in the tray.

The product has been tested using different work platforms. For your information, the commercial names of the work platforms used are specified, which does not imply an exclusive use of such trademarks but the work protocol. The product works correctly as long as the platform used ensures the temperature program raised for the product.

→ See the Strip developing diagram in the Annex of the document.

3.- Strip interpretation

Strip interpretation can be done visually, using the evaluation chart included in the kit (see working procedures for each instrument) or automatically, with BLOTRix R2 or BLOTRix S1 scanners (BioSciTec).

On request, OPERON, S.A. will provide the user with the instructions for the automatic strip interpretation with the scanners (DO-09053004 "Instructions of use BLOTRix R2_S1"). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>


Important: results of strips interpretation with the scanners must always be visually verified; very faint bands could be not detected by the instruments.

RESULTS

For each gene, three potential results may be obtained:

1. **Normal homozygous:** a single blue band appears in the position associated to gene normal sequence. Both chromosomes contain the normal sequence of the gene (N/N).
2. **Normal/mutated heterozygous:** two blue bands appear in the positions associated with gene normal and mutated sequences. One chromosome contains a gene with the normal sequence, and the other chromosome contains a gene with the mutated sequence (N/M).
3. **Mutated homozygous:** a single blue band appears in the position associated to gene mutated sequence. Both chromosomes contain the mutated sequence of the gene (M/M).

Example of genotyping using ThromboStripquality control

QUALITY CONTROL

The development control line should always appear above the upper black line. Its absence will be indicative of issues encountered during the product hybridization/development stage.

The strip corresponding to the amplification negative control will show only the developing control band. Regarding to DNA samples, everyone must present at least one blue band for each gene (except in the case of a potential gene deletion), and two at the most. The absence of bands associated to the genes is indicative of a PCR failure (check DNA quality, verify PCR in a 3% agarose gel) or a hybridization failure (check the hybridization temperature).

PERFORMANCE

The use of this test is restricted to professional users that are familiar with the methods used in Molecular Biology.

The test has been evaluated with 156 DNA samples and the results have been 100 % consistent with those obtained by other genotyping standard methods (restriction assay, sequencing).

POTENTIAL ISSUES

1. No band appears on the strip, including the control band.

- Conjugate and/or developer have not been brought to room temperature.
- Conjugate and/or developer have not been added, or the amount added is insufficient.

2. Only the control band appears.

- PCR amplification failure (verify in agarose gel).
- PCR has not been added, or the amount added is insufficient, at the hybridization step.
- Improper Hybridization Buffer/Wash Buffer 1 temperature (above the indicated).
- Improper incubator temperature (above the indicated).

3. Strong background color on the strip following developing.

- Washing steps were not properly performed (inadequate time, insufficient buffer amount, cold buffers).

4. Irregular strip dyeing.

- Improper agitation (check the plate thermoshaker program).
- Strips were not completely immersed during the incubations.

5. Unexpected results.

- Improper incubation temperature and/or buffer/reagent temperature.
- PCR contamination (check using a negative control).
- Cross-contamination of adjacent channels due to fluid overflow from one channel to the other when adding hybridization buffer.

The following table shows in grey the appropriate temperature working ranges (hybridization/developing; only tested with PST-60HL thermoshaker) for each probe of ThromboStrip (good hybridization intensity, no cross reaction):

	39 °C	40 °C	41 °C	42 °C	43 °C	44 °C	45 °C
PTB N							
PTB M							
FVL N							
FVL M							
MTHFR N							
MTHFR M							

SENSITIVITY

The sensitivity of ThromboStrip was evaluated with several lots of product. The minimum DNA quantity per PCR that can be right genotyped is 5-10 ng (bands with good intensity). DNA quantities lower than 5 ng cause weak bands, in some cases difficult to see.

SPECIFICITY

ThromboStrip specificity was evaluated with 56 DNA samples, of a known genotype, with three different lots of product. In all of them, the results concordance with the previous genotype was of 100 %.

ThromboStrip specificity was evaluated with 100 DNA samples, of a known genotype, with one lot of product. In all of them, the results concordance with the previous genotype was of 100 %.

HOOK EFFECT

ThromboStrip results were evaluated with increased quantities of DNA in the PCR reaction (up to 5000 ng). Neither PCR inhibition nor intensity bands reduction nor crossreactions were observed. The results were right with all the DNA quantities tested.

INTRA-ASSAY PRECISION

5 replica of 3 different DNA samples were analyzed with three lots of ThromboStrip. The study was made by one person and the replica results were always the same, showing a high intra-assay precision result for ThromboStrip.

INTER-DAY PRECISION

Three different DNA samples were analyzed, in duplicate, with three lots of ThromboStrip for 5 consecutive days. The study was made by one person and the results for each one of the DNA samples were the same everyday, showing a high inter-day precision result for ThromboStrip.

INTER-LABORATORY PRECISION

Three DNA samples were analyzed, in duplicate, with the same lot of ThromboStrip by three people. The study was made in one day and the results were the same for all the DNA samples, showing a high inter-laboratory precision for ThromboStrip kit.

INTER-LOT PRECISION

Three different DNA samples were analyzed, in triplicate, with three lots of ThromboStrip. The study was made in one day by one person and the results were allways the same, showing a high inter-lot precision for ThromboStrip kit.

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1. Bertina R.M. et al. "Mutation in blood coagulation factor V associated with resistance to activated protein C". Nature (1994). Vol 369: 64-67.
2. Poort S.R. et al. "A common genetic variation in the 3-prime-untranslated region of the prothrombin gene is associated with elevated plasma prothrombin levels and an increase in venous thrombosis". Blood (1996). Vol 88: 3698-3703.
3. Frosst P. et al. "A candidate genetic risk factor for vascular disease: a common mutation in methylenetetrahydrofolate reductase". Nature Genet (1995). Vol 10: 111-113.

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OPERON, S.A. - Camino del Plano, 19 - E-50410 Cuarte de Huerva - (Zaragoza) – ESPAÑA