High PapillomaStrip

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HIGH PAPILOMASTRIP

Test for the detection of 19 medium-high risk anogenital HPVs

INTENDED USE

The High PapillomaStrip test is a test based on the reverse blot technique that allows qualitative detection in DNA samples from cervico-uterine smears or biopsies, of 19 human *Papillomavirus* subtypes of medium-high risk: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 69, 73 y 82 (MM4 and IS39).

Take also under consideration other parameters such as symptoms or clinical history when giving the final diagnosis.

PRINCIPLE

High PapillomaStrip kit is based on the principle of reverse hybridization and allows for the detection and identification of 19 medium-high risk genital HPVs in DNA samples from cervico-uterine smears or biopsies.

Human *Papillomavirus* (HPV) is responsible for a variety of skin and mucosa disorders. The clinical manifestations may vary according to the anatomic location of the lesion and the type of HPV involved. To date, over 100 different HPV genotypes have been described, of which more than 40 infect the anogenital tract. Of these, approximately 1/3 is associated with cervical cancer and anal neoplasm. The anogenital HPVs are generally divided into two categories: those with a low oncogenic potential risk (Low-Risk Group) and those with a medium-high oncogenic potential risk (Medium-High-Risk Group). The high-risk HPVs are generally associated with high-grade precancerous lesions and invasive cancer, while low-risk HPVs are frequently found in asymptomatic or benign conditions such as genital warts (1, 2).

The most common types of HPV are:

- Medium-high risk: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 69, 73 and 82.
- Low risk: 6, 11, 34, 40, 42, 43, 44.

The HPV virus contains a single molecule of circular double stranded DNA of approximately 8 Kb. Three regions can be distinguished: a region regulating the transcription and replication processes (URR), a gene region of "early" expression encoding proteins involved in transcription, transformation and replication (E1-E2-E5-E6-E7) and a gene region of "late" expression encoding structural proteins (L1-L2-E4).

Currently, the diagnosis of HPV infection is based on viral DNA detection. The presence of conserved regions in different portions of viral DNA has allowed for the development of consensus primers, such as the systems MY09-MY11 (L1 region), PGMY09-PGMY11 (L1 region), Gp5+-Gp6+ (L1 region), SPF10 (L1 region) and LCR-E7 (E7 region), that enable the detection of a wide spectrum of different HPV genotypes. However, the differences in the oncogenic potential of HPVs make it necessary to identify the type of HPV present in the sample, for which different techniques are used: restriction analysis, sequencing, ELISA, reverse hybridization or microarrays (3, 4, 5, 6, 7, 8).

The procedure of High PapillomaStrip kit consists in three steps: a) DNA extraction, b) amplification by PCR and c) hybridization/developing.

a) DNA extraction

Reagents not included in the kit.

Any standard DNA purification protocol can be used.

b) Amplification by PCR

High PapillomaStrip kit includes all reagents required for the amplification of regions E6-E7 from the 19 HPVs detected by this kit. Since the expression of the E6-E7 genes is directly associated with carcinogenesis, this region has been selected as the detection of these targets may provide significant biological and pathological information. Furthermore, some studies show the existence of a gene deletion process in region L1. This deletion occurs when DNA from high-risk HPV integrates into the epithelial regions, and it has been described that it occurs in approximately 30% of the cases of cervical cancer samples positive for E6 region (9, 10).


Specific primers are used for each of the 19 HPVs in order to avoid problems that may occur with consensus primers for preferred amplification of any particular subtype.

During PCR, the amplification of a control gen (β -globin) also occurs, acting as an indicator of the absence of PCR inhibitors in the amplified DNA sample.

c) Hybridization and developing

At this stage there is a specific binding of DNA fragments amplified during the PCR reaction to a series of probes deposited on nylon membranes.

Probes for each HPV detected are covalently linked to the membrane, along with a probe for the amplification control (β -globin, BGB), a probe for the developing control and two lines, one black and one red, to control the positioning of the strip and as an aid in the interpretation of results (see attached diagram).

The detection of the fragments hybridized to different probes is performed using a conjugate (streptavidin-peroxidase) that binds to a biotin label that is added to the DNA fragments amplified during the PCR. Following the addition of a substrate for the peroxidase (TMB) a blue precipitate generates where hybridization has occurred.

As final result of the test, a band pattern is obtained which is interpreted with the aid of a control strip.

MATERIALS PROVIDED WITH THE KIT

HIGH PAPILOMASTRIP			16-test Kit
	Membranes	STRIPS	16 tests
	8-channel trays	PL	2
PCR Reagents REAG PCR	PCR Mix	REAG PCR MIX	0.8 ml
	Primers	REAG PRIMER	0.11 ml
	Taq	REAG TAQ HS	0.050 ml
	Denaturant	SOLN DN	1 ml
	Hybridization Buffer	BUF HYB	60 ml
	Wash Buffer 1	BUF WASH 1	100 ml
	Conjugate	CONJ HRP	60 ml
	Wash Buffer 2	BUF WASH 2	130 ml
	Substrate	SUBS TMB	30 ml
	Instructions of use		1
	Interpretation chart		1

MATERIALS NOT INCLUDED IN THE KIT

The following additional material is required when using the kit:

1. Microtubes for PCR
2. Micropipettes and micropipette tips (sterile or UV-irradiated and ideally with a filter)
3. Tweezers and a pencil (optional)
4. Chronometer
5. Thermometer

REQUIRED EQUIPMENTS FOR KIT DEVELOPMENT

1. Thermoshakers for plates
The kit has been validated with the following thermoshakers: PST-60HL (Biosan S.L.), Profiblot T30 y T48 (Tecan), Autoblot 3000H (medTEC), TENDIGO (Fujirebio), Dynablot Heat (Dynex).
2. Thermocycler
The following thermocyclers have been successfully used with the kit: PTC-100 (MJ Research), MJ Mini Gradient Thermal Cycler (BioRad), Mastercycler Personal (Eppendorf), S1000 Thermal Cycler (BioRad).
3. BLOTriX S1 or BLOTriX R2 scanner (BioSciTec, optional).

PRECAUTIONS

1. Only use the reagents *in vitro*.
2. Strictly follow the instructions provided. Modifying any of the prescribed **steps or temperatures** may severely affect the test results.
3. All the materials and reagents used must be free of DNAses. The use of filtered pipette tips is recommended for PCR preparation to prevent aerosol contamination. Ensure that all the usual precautions are taken in the preparation of the amplifications. Use autoclaved material for the hybridization/development process.
4. Store the kit components as indicated in the instructions.
5. Do not exchange components from kits with different lot numbers.
6. Do not use kit components after the expiration dates.

7. Strips are for single-use.
8. Supplied trays are for single-use.
9. Each channel of the supplied trays is for just one strip.
10. If the package is broken, the product can still be used providing none of the components have been damaged.
11. The used product should be discarded in compliance with current legislation.
12. Patient samples must always be treated as potentially infectious. Environmental and safety standards must be adhered to.
13. Once it is used, it is advisable to manipulate the strip taking into account the same considerations as for the sample. The strips have to be manipulated with the proper precautions and should be managed as biohazardous material.
14. The denaturation solution contains < 2% NaOH and is an irritant to both eyes and skin (H314 and P280, P305, P351, P338, P310).
15. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking and component lots.

STORAGE

Store all reagents at 2-8 °C. Precipitate formation may occur in some solutions (hybridization buffer, wash buffer 1 and 2, conjugate) due to storage at 2-8 °C. These precipitates are redissolved when brought to room temperature (conjugate) or heated to 42 °C hybridization buffer and wash buffer 1 and 2).

Expiry date for all reagents is printed on the respective label.

SAMPLES

The test has been designed and validated for its use with DNA derived from cervico-uterine smears. It has also shown to perform correctly with samples from penis swabs.

For sampling use dry and sterile cotton swabs or brushes. Make sure to take sufficient sample amount, but without causing bleeding from the lesion. Put the swab into its tube, do not use any preserving medium. Store at 2-8 °C if it is to be processed within one week or at -20 °C for a later processing.

The quality and concentration of DNA extracted from cervico-uterine specimens is crucial for the proper functioning of the test. For this reason, the kit includes a control that detects the presence of PCR inhibitors in DNA (amplification of the β -globin). Additionally, we recommend measuring the concentration of DNA samples to be used, since the use of DNA samples with a concentration below 10 ng/ μ l, although showing control gene amplification, could reduce the sensitivity of the test.

Some of the DNA extraction procedures that have been evaluated with satisfactory results with High PapillomaStrip kit are:

- a) DNA extraction kit from OPERON S.A.: Ref 3.149.020.64.000, DNA extraction kit.
- b) DNA extraction with silica resins (QIAamp DNA Mini Kit; Qiagen).
- c) Extraction protocol for microarray Clart HPV2 from Genomica S.L.
- d) DNA extraction with magnetic particles.

Store the DNA samples at 2-8 °C if they are going to be used shortly after, or at -20 °C for longer periods of storage.

HIGH PAPILOMASTRIP PROCEDURE

1.- Polymerase chain reaction

PCR preparation

Important: before opening the vials with the PCR reagents centrifuge them briefly. This will ensure that all the contents will be at the bottom of the tube.

- Place the tubes required according to the number of DNAs to be amplified.
- Add to each PCR tube: 38 μ l of PCR premix + 5 μ l of primers + 2 μ l of Taq + 5 μ l of DNA. Mix. Whenever possible, keep all reagents at 2-8 °C during preparation. If several DNAs are to be amplified, it is recommended to prepare a common mixture containing all reagents, so 45 μ l mixture + 5 μ l DNA are finally added. For example:

Nr of PCRs	PCR Premix	Primers	Taq
1	38 μ l	5 μ l	2 μ l
3	133 μ l	17.5 μ l	7 μ l
5	209 μ l	27.5 μ l	11 μ l
8	342 μ l	45 μ l	18 μ l

* Always prepare an excess of the mixtures containing all PCR reagents, to compensate volume losses during pipetting process.

Due to requirements of "good laboratory practice", the user should also include a negative control (e.g. water or TE buffer) to exclude contamination and, if required, a positive control (not included in the kit).

Amplification

Regard the importance of getting a suitable concentration of the DNA to get a proper result.

Place the tubes into the cycler device (if required, add 1 drop of Nujol oil on top of every PCR reaction) and amplify DNA by the following program:

94 °C , 5 min
 40 cycles of: 94 °C, 1 min
 58 °C, 1 min
 72 °C, 1 min
 72 °C, 5 min
 4 °C

Once the PCR is finished, continue with the developing of the strip. If the samples are not to be tested immediately, they must be stored at 2-8 °C for no more than 24 hrs. For longer storage, freeze at -20 °C.

→ See the PCR preparation diagram in the Annex of the document.

2.- Strip developing

On request, Operon will provide the user with the working procedure for the specific instrument to be used for processing the strips (PST-60HL, Profiblot, Autoblot 3000H, TENDIGO or Dynablot Heat). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>

In the denaturing step, **12.5 μ l of denaturation buffer + 12.5 μ l of PCR** will be dispensed in the tray.

The product has been tested using different work platforms. For your information, the commercial names of the work platforms used are specified, which does not imply an exclusive use of such trademarks but the work protocol. The product works correctly as long as the platform used ensures the temperature program raised for the product.

→ See the Strip developing diagram in the Annex of the document.

3.- Strip interpretation

Strip interpretation can be done visually, using the evaluation chart included in the kit (see working procedures for each instrument) or automatically, with BLOTRix R2 or BLOTRix S1scanners (BioSciTec).

On request, OPERON, S.A. will provide the user with the instructions for the automatic strip interpretation with the scanners (DO-09053004 "Instructions of use BLOTRix R2_S1"). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>

Important: results of strips interpretation with the scanners must always be visually verified; very faint bands could be not detected by the instruments.


RESULTS

High PapillomaStrip kit allows for the detection and genotyping of 19 of the most common medium-high risk anogenital HPVs: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 69, 73 and 82. The bands obtained on each strip are interpreted with the aid of a chart.

A DNA sample may render the following results:

- The control band only (β -globin): none of the 19 HPVs detected by this kit is present or could be amplified.
- The control band and one band associated to a certain type of HPV: sample positive for that particular HPV.
- The control band and several bands associated to different HPVs: positive sample with multiple infection.
- One or several bands associated to different HPVs: positive sample with infection by one or several HPVs in which the β -globin gene does not amplify. This may happen, for example, when the sample has a high number of copies of one or more viruses, as the PCR conditions have been designed to favor the amplification of HPV DNA versus that of the internal control.

Examples of genotyping using High PapillomaStrip

QUALITY CONTROL

The development control line should always appear above the upper black line. Its absence will be indicative of issues encountered during the product hybridization/development stage.

The strip corresponding to the amplification negative control will show only the developing control band. For negative samples, the presence of the band associated to the internal control (β -globin) is

required to confirm that it is a true negative. The absence of this band in a negative sample is indicative of the absence of amplification during the PCR, which can be caused by lack of DNA (measure concentration), by degradation or by the presence of PCR inhibitors in the same (hemoglobin , traces of salts, etc.).

FEATURES / EXTERNAL AND INTERNAL EVALUATIONS

The use of this test is restricted to professional users that are familiar with the methods used in Molecular Biology.

The test has been evaluated with several DNA samples from cervico-uterine smears, obtained by different protocols, with satisfactory results. The test has been successfully used for the analysis of the presence of HPV in penis swabs.

Evaluation nº 1

The sensitivity and specificity of High PapillomaStrip were compared with the ones of HC2 HPV DNA test from Digene/Qiagen. For this purpose, a total of 50 samples of cervico-vaginal smears were analyzed by both methods and the results were verified by amplification with individual primers specific to each one of the HPVs detected, followed by agarose gel analysis and hybridization with specific probes.

Cervico-vaginal smears were processed according to the protocol of HC2 HPV DNA test (denaturation by addition of NaOH). Once performed the Digene test, we proceeded to the purification of DNA from these samples by using silica columns (QuiAprep Spin MiniPrep kit, Qiagen). These DNA samples were directly used for the analysis using three different High PapillomaStrip lots.

The results of sensitivity and specificity concordance are shown in the table below:

		+ (42)	- (8)	Total
High PapillomaStrip	+	41	0	41
	-	1	8	9
	Total	42	8	50
Sensitivity Concordance = 97.6 % Specificity Concordance > 99 % Concordance Percentage = 98 %				

		+ (42)	- (8)	Total
Digene	+	31	0	31
	-	11	8	19
	Total	42	8	50
Sensitivity Concordance = 73.8 % Specificity Concordance > 99 % Concordance Percentage = 78 %				

Evaluation nº 2

The sensitivity and specificity of High PapillomaStrip were compared with the Human *Papillomavirus* CLART test from Genomica (Zeltia Group). The study was performed at the Centro de Análisis Genéticos C.A.G.T. (Zaragoza). A total of 99 samples of cervico-vaginal smears were analyzed by both methods and the results were verified by amplification with individual primers specific to each one of the HPVs detected, followed by agarose gel analysis and hybridization with specific probes.

DNA was extracted from the cervico-vaginal smears according to the protocol of the Human Papillomavirus CLART kit. These DNA samples were used both for the analysis using the microarray and High PapillomaStrip.

The results on sensitivity and specificity concordance are shown in the table below:

		+ (59)	- (40)	Total
High PapillomaStrip	+	57	0	57
	-	2	40	42
	Total	59	40	99
Sensitivity Concordance = 96,6 % Specificity Concordance > 99 % Percentage of Concordance = 98 %				

		+ (59)	- (40)	Total
CLART	+	53	0	53
	-	6	40	46
	Total	59	40	99
Sensitivity Concordance = 89,8 % Specificity Concordance > 99 % Percentage of Concordance = 94 %				

The percentage of detection and results on sensitivity concordance of each HPV subtype within the set of samples tested was:

HPV	Detection percentage in sample		Sensitivity concordance		HPV	Detection percentage in sample		Sensitivity concordance	
	High PapillomaS trip	Microarray CLART HPV	High Papilloma Strip	Microarray CLART HPV		High Papilloma Strip	Microarray CLART HPV	High Papilloma Strip	Microarray CLART HPV
16	16,1	14,1	94,1	82,3	53	5	6,1	71,4	85,7
18	3	3	100	100	56	12,1	6,1	100	50
26	0	0	--	--	58	7,1	9,1	77,8	100
31	5	5	100	100	59	7,1	5	100	71,4
33	4	4	100	100	66	13,1	13,1	92,9	92,9
35	3	2	100	66,7	68	7,1	3	100	42,9
39	2	0	100	0*	69	0	0	--	--
45	0	0	--	--	73	4	1	100	25
51	6,1	7,1	85,7	100	82	4	4	100	100
52	9,1	7,1	100	77,8					

*Result obtained with only two positive samples.

POTENTIAL ISSUES

1. No bands appear on the strip, including the control band.

- The conjugate and/or developer were not balanced at room temperature.
- The conjugate and/or developer were not added, or were added in too small quantity.

2. Only the control band appears

- PCR amplification has failed (check with an agarose gel).
- The PCR was not added, or was added in too small quantity in the hybridization step.
- Incorrect hybridization or wash 1 buffers temperature (higher than instructed).
- Incorrect incubator temperature (higher than instructed).

3. The strip has a strong background colour following development

- The washing steps were not performed effectively (wrong timing, buffer quantity too small, cold buffers).

4. Non-homogenous strip staining

- Inadequate shaking (please review the plate thermo-shaker programme).
- The strips were not entirely submerged during the incubations.

5. Unexpected results

- Inadequate incubation and/or buffer/reagent temperatures.
- PCR contamination (verify using a negative control).
- Contamination of adjacent channels owing to the transfer of liquid from one channel to another when adding the hybridization buffer.

The following table displays in gray the working temperature ranges (hybridization/developing; test performed using the PST-60HL thermoshaker) with which adequate results are obtained for each band (good hybridization intensity, no cross-reaction):

	39 °C	40 °C	41 °C	42 °C	43 °C	44 °C	45 °C
SBG							
S16							
S18							
S26							
S31							
S33							
S35							
S39							
S45							
S51							

	39 °C	40 °C	41 °C	42 °C	43 °C	44 °C	45 °C
S52							
S53							
S56							
S58							
S59							
S66							
S68							
S69							
S73							
S82							

SENSITIVITY

The detection sensitivity of High PapillomaStrip for each of the HPVs detected by the kit was determined.

For this purpose, plasmids of regions E6-E7 from the 19 HPVs were prepared in the pGemT vector. The concentration of each plasmid was quantified by spectrophotometry at 260 nm and the High PapillomaStrip detection sensitivity for each HPV was determined by the analysis, using three different lots of product, of 1/10 dilutions of the different plasmids in a negative control DNA.

Sensitivity results obtained are shown in the following table:

HPV Subtype	Sensitivity	
	Concentration (pg/ µl)	pg per PCR
HPV16	0,0002	0,001
HPV18	0,0002	0,001
HPV26	0,00002	0,0001
HPV31	0,00002	0,0001
HPV33	0,00002	0,0001
HPV35	0,00002	0,0001
HPV39	0,0002	0,001
HPV45	0,00002	0,0001
HPV51	0,0002	0,001
HPV52	0,00002	0,0001

HPV Subtype	Sensitivity	
	Concentration (pg/ µl)	pg per PCR
HPV53	0,00002	0,0001
HPV56	0,00002	0,0001
HPV58	0,0002	0,001
HPV59	0,0002	0,001
HPV66	0,0002	0,001
HPV68	0,0002	0,001
HPV69	0,0002	0,001
HPV73	0,00002	0,0001
HPV82	0,002	0,01

SPECIFICITY

The absence of reaction was checked with the following pathogens:

Escherichia coli	HPV6	HPV74
Neisseria gonorrhoeae tipo I	HPV11	HPV81
Neisseria gonorrhoeae tipo II/III	HPV40	HPV83
Herpes simplex 1	HPV42	HPV84
Herpes simplex 2	HPV43	HPV91

Haemophilus ducreyi	HPV44	
Trichomonas vaginalis	HPV54	
Streptococcus agalactiae	HPV61	
Staphylococcus aureus mecA (-)	HPV62	
Staphylococcus aureus mecA (+)	HPV67	
Cándida albicans	HPV70	
Ureaplasma urealyticum	HPV71	
Mycoplasma genitalium	HPV72	

HOOK EFFECT

The performance of High PapillomaStrip kit in the presence of increasing amounts of plasmidic DNA (up to 100 ng) was evaluated. No PCR inhibition effect, or reduction in the intensity of the bands obtained, or problems arising from cross-reaction were observed. Therefore, there is no hook effect.

INTRA-ASSAY PRECISION

The intra-assay precision of High PapillomaStrip has been studied examining, first, 5 sensitivity curves (1/10 dilutions of a control DNA consisting of a mixture of plasmids of the 19 HPVs detected by the kit) and on the other side, five replicates of five DNA samples obtained from cervico-vaginal smears.

Overall, High PapillomaStrip kit has shown a good intra-assay precision. Both the test performed with dilution curves (where all detected HPVs are amplified simultaneously) and the test using actual samples, show a very good reproducibility for samples of medium, normal or strong intensity and a slight variation in samples near the limit of detection of the different HPVs.

INTER-DAY PRECISION

The inter-day precision of High PapillomaStrip has been studied examining 5 sensitivity curves (1/10 dilutions of a control DNA consisting of a mixture of plasmids of the 19 HPVs detected by the kit). Each sensitivity curve was constructed by the same operator, using the same product lot and over five consecutive days.

Overall, High PapillomaStrip kit shows a good reproducibility varying the day. As in the intra-assay precision, a very good reproducibility for samples of medium, normal or strong intensity and a slight variability for the samples near the limit of detection of the different HPVs can be observed.

INTER-LAB PRECISION

The inter-lab precision has been studied examining 3 sensitivity curves (1/10 dilutions of a control DNA consisting of a mixture of plasmids of the 19 HPVs detected by the kit). Each curve was constructed by a different operator, on the same day and using the same product lot.

Overall, the High PapillomaStrip kit shows a good reproducibility varying the operator and the conditions. The largest differences are found, as for the rest of the tests performed with these samples (mixtures of all detected HPVs), in the dilutions near the limit of sensitivity of each HPV.

INTER-LOT PRECISION

The inter-lot precision of High PapillomaStrip kit has been studied by the analysis using three different product lots of a sensitivity curve (1/10 dilutions of a control DNA consisting of a mixture of plasmids of the 19 HPVs detected by the kit) and 19 plasmid samples from each of the HPVs detected by this kit.

Overall, High PapillomaStrip kit has shown a good inter-lot precision. The results obtained for the 19 plasmid samples were identical in all cases, the results of the sensitivity curves, as in the rest of the precision tests, have demonstrated a very good reproducibility for samples of medium, normal or strong intensity and a light variability for the samples near the limit of detection of the different HPVs.

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1. Bosch, F.X. et al. "Prevalence of human papillomavirus in cervical cancer: a worldwide perspective. International biological study on cervical cancer (IBSCC) Study Group". J. Natl. Cancer Inst. (1995). 87: 796-802.
2. Bosch, F.X. et al. "Papillomavirus research update: highlights of the Barcelona HPV 2000 international papillomavirus conference". J. Clin. Pathol. (2001). 54: 163-175.
3. Gravitt P.E. et al. "Genotyping of 27 human papillomavirus types by using L1 consensus PCR products by a single-hybridization, reverse line blot detection method". J. Clin. Microbiol. (1998). Vol 36, N° 10: 3020-3027.
4. Coutlèe F. et al. "Enhanced detection and typing of human papillomavirus (HPV) DNA in anogenital samples with PGMV primers and the linear array HPV genotyping test". J. Clin. Microbiol. (2006). Vol 44, N° 6: 1998-2006.
5. Jacobs M.V. et al. "Group-specific differentiation between high- and low-risk human papillomavirus genotypes by general primer-mediated PCR and two cocktails of oligonucleotide probes". J. Clin. Microbiol. Vol 33, N° 4: 901-905.
6. Kleter B. et al. "Development and clinical evaluation of a highly sensitive PCR-reverse hybridization line probe assay for detection and identification of anogenital human papillomavirus". J. Clin. Microbiol. (1999). Vol 37, N° 8: 2508-2517.
7. Gheitt T. "Development of a sensitive and specific assay combining multiplex PCR and DNA microarray primer extension to detect high-risk mucosal human papillomavirus types". J. Clin. Microbiol. (2006). Vol 44, N° 6: 2025-2031.
8. Sandri M.T. "Comparison of the digene HC2 assay and the Roche AMPLICOR human papillomavirus (HPV) test for detection of high-risk HPV genotypes in cervical samples". J. Clin. Microbiol. (2006). Vol 44, N° 6: 2141-2146.
9. Han, J. et al. "Simultaneous amplification and identification of 25 human papillomavirus types with Templex Technology". (2006). J. Clin. Microbiol. Vol 44, N°11: 4157-4162.
10. Sotlar K. et al. "Detection and typing of human papillomavirus by E6 nested multiplex PCR". (2004). J. Clin. Microbiol. Vol 42, N° 7: 3176-3184.

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OPERON, S.A. - Camino del Plano, 19 - E-50410 Cuarte de Huerva - (Zaragoza) – SPAIN