STD Panel Strip

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STD Panel Strip

Test for the detection of 10 pathogens associated to sexually transmitted diseases

INTENDED USE

The STD Panel Strip test is a test based on the reverse blot technique that allows the qualitative detection of 10 pathogens associated to sexually transmitted diseases: *Chlamydia trachomatis* (discriminating variants L1, L2 and L3, that cause lymphogranuloma venereum), *Neisseria gonorrhoeae*, *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Ureaplasma parvum*, *Ureaplasma urealyticum*, *Mycoplasma hominis*, *Herpes simplex 1*, *Herpes simplex 2*, and *Treponema pallidum*.

Through comparative analysis of described sequences using BLAST and BioEdit software the detection of all different serotypes/genetic sequences described for all the pathogens has been verified. So, all *Chlamydia trachomatis* serotypes (A, B, B_α, C, D, D_α, E, F, G, H, I, I_α, J, J_α, K, L₁, L₂ y L₃) and *Ureaplasma* spp serovars (1-14) are detected.

The product may be used for any indication that involves the identification of the previously mentioned pathogens: initial diagnosis, screening, monitoring.

Take also under consideration other parameters such as symptoms or clinical history when giving the final diagnosis.

PRINCIPLE

The STD Panel Strip kit is based on the principle of reverse hybridization, and allows the detection and identification of 10 pathogens associated to sexually transmitted diseases in DNA samples.

At present, sexually transmitted diseases (STDs) make up the most common group of infectious diseases that must be reported in the majority of the world's countries. Their incidence is elevated, and more than 20 pathogenic micro-organisms are known to be transmitted by sexual contact, with *Chlamydia trachomatis*, *Neisseria gonorrhoeae*, *Trichomonas vaginalis* and *Mycoplasma* (*genitalium* and *hominis*) being the most common.

The majority of STDs have treatments and are curable. They affect both men and women, with 2/3 of cases in those under 25 years of age. They can occur without any symptoms appearing, especially in women, but the disease is not cured unless appropriate treatment is received. If untreated, they can cause serious health problems, such as infertility, blindness, mental disorders, physical birth defects, increased odds of developing cancer, heart disease, and even death. Depending on the disease, the infection can be transmitted through any type of activity involving the sexual organs; it can also be transmitted by contact with blood during sexual activity.

The current permanent trend of increasing sexually transmitted diseases, especially among young people under 25 years, makes their early and accurate diagnosis of great interest.

The procedure of kit STD Panel Strip consists in three steps: a) DNA extraction, b) amplification by PCR and c) hybridization/developing.

a) DNA extraction

Reagents not included in the kit.
See "SAMPLES" section.

b) PCR amplification

The STD Panel Strip includes the reagents needed to perform specific genetic sequence amplification of the pathogens detected with the kit. The genes used in each case are presented in the following table:


| Pathogen | Gen / Locus | Pathogen | Gen / Locus |
|-----------------------|-----------------|---------------------------------------|------------------|
| <i>C. trachomatis</i> | Omp | <i>Herpes virus</i> | DNA polimerase |
| | Cryptic plasmid | <i>C. trachomatis</i> (L1, L2, L3) | pmpH |
| <i>T. vaginalis</i> | TRIDNATARP | <i>T. pallidum</i> | DNA polimerase I |
| <i>M. genitalium</i> | MgPA 190 | Control gen B | β -globin |
| <i>Ureaplasmas</i> | Urease (Ure B) | | |
| <i>M. hominis</i> | gap | | |
| <i>N. gonorrhoeae</i> | Opa | | |
| | Por A | | |
| Control gen A | β -globin | | |

During PCR, two regions of the human control gene (β -globin) are also amplified.

c) Hybridization and development

Specific binding of PCR-amplified DNA fragments takes place during this step, using a series of probes that are covalently bound to a nylon membrane.

The membranes are covalently bound to probes specific to each one of the detected microorganisms, along with probes for the two regions of the control gene (BGB and BGBO1), a developing control line and two lines, one black and one red, for the control of the strip's position to aid in the interpretation of results (see attached drawing).

The detection of the fragments that hybridize to the different probes is carried out using a streptavidin-peroxidase conjugate that binds to a biotin marker added to the PCR-amplified DNA fragments. Following the addition of a peroxidase substrate (TMB), a blue precipitate appears in the locations where hybridization has taken place. The end result of the test is a pattern of bands that can be interpreted with the help of a control strip.

KIT CONTENTS

| STD Panel Strip | | | Kit 16 tests | Kit 48 tests |
|--------------------------|----------------------|--------------|------------------|------------------|
| | Membranes | STRIPS | 16 strips | 48 strips |
| | 8 channel tray | PL | 2 | -- |
| PCR reagents REAG PCR | PCR Mix | REAG PCR MIX | 2 x 0,80 ml | 3 x 1,40 ml |
| | Primers | REAG PRIMER | 0,110 + 0,110 ml | 0,275 + 0,275 ml |
| | Taq | REAG TAQ HS | 0,050 ml | 0,120 ml |
| | Denaturing | SOLN DN | 1 ml | 2 x 1 ml |
| | Hybridization buffer | BUF HYB | 60 ml | 125 ml |
| | Wash buffer 1 | BUF WASH 1 | 100 ml | 250 ml |
| | Conjugate | CONJ HRP | 60 ml | 125 ml |
| | Wash buffer 2 | BUF WASH 2 | 130 ml | 350 ml |
| | Substrate | SUBS TMB | 30 ml | 65 ml |
| | Instructions | | 1 | 1 |
| | Interpretation chart | | 1 | 1 |

MATERIALS NOT INCLUDED IN THE KIT

The following additional material is required when using the kit:

1. Microtubes for PCR
2. Micropipettes and micropipette tips (sterile or UV-irradiated and ideally with a filter)
3. Tweezers and a pencil (optional)
4. Chronometer
5. Thermometer

REQUIRED EQUIPMENTS FOR KIT DEVELOPMENT

1. Thermocycler
The following thermocyclers have been successfully used with the kit: 2720 Thermal Cycler (Applied Biosystems), MJ Mini Gradient Thermal Cycler (BioRad), Mastercycler Personal (Eppendorf), S1000 Thermal Cycler (BioRad).
2. Thermoshakers
The kit has been validated with the following thermoshakers: PST-60HL (Biosan S.L.), Profiblot T30 and T48 (Tecan), Autoblots 3000H (medTEC), TENDIGO (Fujirebio), Dynablot Heat (Dynex).
3. BLOTrix S1 or BLOTrix R2 scanner (BioSciTec, optional).

PRECAUTIONS

1. Only use the reagents *in vitro*.
2. Strictly follow the instructions provided. Modifying any of the prescribed **steps or temperatures** may severely affect the test results.

3. All the materials and reagents used must be free of DNAses. The use of filtered pipette tips is recommended for PCR preparation to prevent aerosol contamination. Ensure that all the usual precautions are taken in the preparation of the amplifications. Use autoclaved material for the hybridization/development process.
4. Store the kit components as indicated in the instructions.
5. Do not exchange components from kits with different lot numbers.
6. Do not use kit components after the expiration dates.
7. Strips are for single-use.
8. Supplied trays are for single-use.
9. Each channel of the supplied trays is for just one strip.
10. If the package is broken, the product can still be used providing none of the components have been damaged.
11. The used product should be discarded in compliance with current legislation.
12. Patient samples must always be treated as potentially infectious. Environmental and safety standards must be adhered to.
13. Once it is used, it is advisable to manipulate the strip taking into account the same considerations as for the sample. The strips have to be manipulated with the proper precautions and should be managed as biohazardous material.
14. The denaturation solution contains < 2% NaOH and is an irritant to both eyes and skin (H314 and P280, P305, P351, P338, P310).
15. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking and component lots.

STORAGE

Store all reagents at 2-8 °C. Precipitates may appear in some solutions (hybridization buffer, wash buffers 1 and 2, conjugate) owing to storage at 2-8 °C. The precipitates are dissolved again when they are brought back to room temperature (conjugate) or heated to 42 °C (hybridization buffer and wash buffers 1 and 2).

The expiry date of all the reagents is printed on the label.

SAMPLES

The test has been designed and validated for use with DNA obtained from smears of various origins (urethral, endocervical, rectal, vaginal), and from urine.

For the smears, use sterile and dry cotton swabs or brushes to take the sample. Ensure enough sample is obtained, but without causing bleeding. Place the swab in a tube in an adequate transport medium (for example, Amies-Stuart medium or Cobas PCR Female Swab Sample Packet by Roche Molecular Systems).

For urine samples, you should ask the patient to collect the first 10-20 ml of urine in a sterile tube, having refrained from urinating for at least two hours before sampling.

The quality and concentration of the DNA extracted is of vital importance for the proper performance of the kit. The presence of inhibitors in the sample or low DNA concentration could reduce the sensitivity of the test.

Some of the DNA extraction procedures that have been evaluated with satisfactory results with STD Panel Strip are:

- a) Manual DNA extraction with silica resins: QIAamp DNA Mini Kit (it is recommended to follow the procedure described for viruses) or MinElute Virus Spin Kit (Qiagen).
In the case of dry swabs, resuspend its content with 200 – 1000 µl of phosphate buffer (PBS) and continue with the extraction.
With QIAamp DNA Mini Kit: we recommend the extraction of 200 µl of smear and elution with 100 µl of TE or elution buffer.

With MinElute Virus Spin Kit: we recommend the extraction of 200 μl of smear and elution with 20-50 μl of TE or elution buffer.

- b) DNA extraction with COBAS 4800 platform (Roche Diagnostics): follow the manufacturer's instructions.
- c) DNA extraction with QiaCube platform (Qiagen): follow the manufacturer's instructions.

Store the DNA samples at 2-8 °C if they are going to be used shortly after, or at -20 °C for longer periods of storage.

INTERFERING SUBSTANCES

Cervical specimens may contain interfering substances that are known to inhibit PCR reactions. These substances include, but are not limited to, the following: whole blood (human), leukocytes and contraceptive and feminine hygiene products, like contraceptive jellies, anti-fungal creams, spermicides, vaginal lubricants, etc).

Most of these interfering substances will be eliminated during DNA extraction. In case of a DNA sample with any inhibitor, it will be detected by the absence of amplification of control genes amplified during the PCR reaction.

STD PANEL STRIP PROCEDURE

1.- Polymerase chain reaction

PCR preparation

Important: before opening the vials with the PCR reagents centrifuge them briefly. This will ensure that all the contents will be at the bottom of the tube.

For each sample two amplification reactions will be run: one of them with primers A (vial with blue cap) and another one with primers B (vial with white cap). PCR Premix and Taq reagents are common to both reactions.

- Prepare the required PCR tubes based on the number of DNAs to be amplified (two per sample).
- Add to each PCR tube: 39 μl of PCR premix + 5 μl of primers (A or B) + 1 μl of Taq + 5 μl of DNA. Mix well. If possible, keep the reagents and the mixtures at 2-8 °C during the preparation. If various DNA samples are to be amplified, preparing a common mixture with all the reagents is recommended to finally add 45 μl of mixture + 5 μl of DNA. For example:

| N° of PCRs | PCR premix | Primers (A or B) | Taq |
|------------|---------------------|--------------------|-------------------|
| 1 | 39 μl | 5 μl | 1 μl |
| 3 | 136,5 μl | 17,5 μl | 3,5 μl |
| 5 | 214,5 μl | 27,5 μl | 5,5 μl |
| 8 | 351 μl | 45 μl | 9 μl |

- *The mixtures containing all PCR reagents should always be prepared in excess to compensate for the loss of volume that takes place during the pipetting process.*

Due to requirements of "good laboratory practice", the user should also include a negative control (e.g. water or TE buffer) to exclude contamination and, if required, a positive control (not included in the kit).

Amplification

Regard the importance of getting a suitable concentration of the DNA to get a proper result.

Place the tubes into the cycler device (if required, add 1 drop of Nujol oil on top of every PCR reaction) and amplify DNA by the following program:

94 °C, 5 min
40 cycles of: 94 °C, 15 sec
57 °C, 1 min 30 sec
72 °C, 45 sec
72 °C, 5 min
4 °C

→ See the PCR preparation diagram in the Annex of the document.

Once the PCR is finished, continue with the developing of the strip. If the samples are not to be tested immediately, they must be stored at 2-8 °C for no more than 24 hrs. For longer storage, freeze at -20 °C.

2.- Strip developing

On request, Operon will provide the user with the working procedure for the specific instrument to be used for processing the strips (PST-60HL, Profiblot, Autoblots 3000H, TENDIGO or Dynablot Heat). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>

For each sample, mix the same volume of PCR A and PCR B (for example, 15 µl + 15 µl) and continue with working procedure for the specific instrument to be used. In the denaturing step, **25 µl of denaturation buffer + 25 µl of the mixture of PCRs** will be dispensed in the tray.

The product has been tested using different work platforms. For your information, the commercial names of the work platforms used are specified, which does not imply an exclusive use of such trademarks but the work protocol. The product works correctly as long as the platform used ensures the temperature program raised for the product.

3.- Strip interpretation

Strip interpretation can be done visually, using the evaluation chart included in the kit (see working procedures for each instrument) or automatically, with BLOTRix R2 or BLOTRix S1 scanners (BioSciTec).

On request, OPERON, S.A. will provide the user with the instructions for the automatic strip interpretation with the scanners (DO-09053004 "Instructions of use BLOTRix R2_S1"). Alternatively, please refer to the following web address:

<https://operon.es/products-protocols/>

Important: results of strips interpretation with the scanners must always be visually verified; very faint bands could be not detected by the instruments.

RESULTS


The STD Panel Strip kit allows the detection and identification of 10 pathogens associated to sexually transmitted diseases: *Chlamydia trachomatis* (CT) and its variants L1, L2 and L3 that cause lymphogranuloma venereum (LGV), *Neisseria gonorrhoeae* (NG), *Mycoplasma genitalium* (MG), *Trichomonas vaginalis* (TV), *Ureaplasma parvum* (UP), *Ureaplasma urealyticum* (UU), *Mycoplasma hominis* (MH), *Herpes simplex 1* (HSV1), *Herpes simplex 2* (HSV2), and *Treponema pallidum* (TP).

A DNA sample may render the following results:

- Only the control bands (BGB and BGBO1): none of the pathogens detected by this kit is present or could be amplified.
- The control bands and one band associated to a certain pathogen: sample positive for that particular pathogen.
- The control bands and several bands associated to different pathogens: positive sample with multiple infection.
- One or several bands associated to different pathogens: positive sample for those pathogens in which the control genes do not amplify. This may happen, for example, when the sample has a high

number of copies of one or more pathogens, as the PCR conditions have been designed to favour the amplification of pathogen DNA versus that of the internal controls. The presence of the control bands is compulsory for a negative sample but optional for a positive sample.

Examples of genotyping STD Panel Strip

Attention: in the case of *Lymphogranuloma venereum* (LGV), as it is a variant of *Chlamydia trachomatis*, the bands associated with *C. Trachomatis* will appear along with those for LGV.

QUALITY CONTROL

The development control line should always appear above the upper black line. Its absence will be indicative of issues encountered during the product hybridization/development stage.

The strip corresponding to the amplification of the negative control will show only the developing control band.

For negative samples, the presence of the band associated to the internal controls (BGB and BGBO1) is required to confirm that it is a true negative. The absence of this band in a negative sample is indicative of the absence of amplification during the PCR, which can be caused by lack of DNA (measure concentration), by degradation or by the presence of PCR inhibitors in the sample (hemoglobin, traces of salts, etc.).

FEATURES / EXTERNAL AND INTERNAL EVALUATIONS

The use of this test is restricted to professional users that are familiar with the methods used in Molecular Biology.

The test has been evaluated with several DNA samples from different origin in several evaluations that are shown below. The first 10 evaluations were performed in the facilities of OPERON, S.A., while the evaluation n° 11 was performed in the premises of the Centro de Análisis Genéticos (CITOGEN, Zaragoza, España). All comercial kits used were CE marked.

Evaluation n° 1

The sensitivity and specificity of detection of *Trichomonas vaginalis* was analysed using 77 DNA samples extracted from fresh cervico-vaginal smears previously evaluated by analysis with light microscopy and confirmation with specific culture medium (Dept. Microbiology, University Hospital of Zaragoza). The extraction of DNA samples was carried out with the DNA Mini Kit (Qiagen).

The concordance of sensitivity and specificity between methods was 100% and 94.4%, respectively. Of the 18 negative samples by microscopy/culture, one gave a positive result with the STD Panel Strip kit. This positive was confirmed by PCR analysis of two alternative genetic regions of the protozoan.

Evaluation n° 2

The sensitivity and specificity of detection of *Chlamydia trachomatis*, serotypes L1, L2 and L3, which cause *Lymphogranuloma venereum*, were analysed. For sensitivity, 10 DNA samples previously evaluated by Real Time PCR (routine procedures at Miguel Servet Hospital in Zaragoza, University Hospital of Donostia in San Sebastián and Valme Hospital in Seville) were used. For specificity, 13 DNA samples extracted with Qiagen's Qiacube system were used and evaluated with the RealCycler THLV kit (Progenie).

The concordance of sensitivity and specificity obtained between methods was 100%.

Evaluation n° 3

The sensitivity and specificity of detection of *Treponema pallidum* were analysed using 13 DNA samples (2 positive and 11 negative) extracted with Qiagen's Qiacube system and evaluated with the RealCycler THLV kit (Progenie).

The concordance of sensitivity and specificity obtained between methods was 100%.

Evaluation n° 4

The sensitivity and specificity of detection of HSV1 were analysed using 24 DNA samples (15 positive and 9 negative) extracted with Qiagen's Qiacube system and analysed with the Fluorotype HSV kit (Hain).

The concordance of sensitivity and specificity obtained between methods was 100%.

Furthermore, the sensitivity and specificity of detection of HSV1 were analysed with 68 samples (28 positive and 40 negative) extracted with the COBAS 4800 platform and analysed with the LightCycler HSV1/2 Qual kit (Roche).

The concordance of sensitivity and specificity obtained between methods was 96.4 and 100%, respectively. A single false negative was found with the STD Panel Strip kit, which corresponded to a sample with a Ct of 40 with the LightCycler HSV1/2 kit (Roche).

Evaluation n° 5

The sensitivity and specificity of detection of HSV2 were analysed using 24 DNA samples (9 positive and 15 negative) extracted with Qiagen's Qiacube system and analysed with the Fluorotype HSV kit (Hain).

The concordance of sensitivity and specificity obtained between methods was 88.9 and 100%, respectively. A positive sample was negative for HSV2 with the STD Panel Strip kit, but positive for *Treponema pallidum* (confirmed with the RealCycler THLV kit by Progenie).

Furthermore, the sensitivity and specificity of detection of HSV1 were analysed with 68 samples (40 positive and 28 negative) extracted with the COBAS 4800 platform and analysed with the LightCycler HSV1/2 Qual kit (Roche).

The concordance of sensitivity and specificity obtained between methods was 97.5 and 100%, respectively. A single false negative was found with the STD Panel Strip kit, which corresponded to a sample with a Ct of 42.5 with the LightCycler HSV1/2 kit (Roche).

Evaluation n° 6

The sensitivity and specificity of detection of *Chlamydia trachomatis* were analysed using 197 DNA samples (55 positive and 142 negative) extracted with Qiagen's Qiacube system and analysed with the Dx CT/NG/MG Assay kit (BioRad).

The concordance of sensitivity and specificity obtained between methods was 100 and 99.3%, respectively. A negative sample tested positive for CT with the STD Panel Strip kit.

Furthermore, the sensitivity of detection of *Chlamydia trachomatis* was analysed with 50 positive samples extracted with the COBAS 4800 platform and analysed with the CT/NG Cobas 4800 kit (Roche).

The concordance of sensitivity obtained between methods was 100%.

Evaluation n° 7

The sensitivity and specificity of detection of *Neisseria gonorrhoeae* were analysed using 193 DNA samples (76 positive and 117 negative) extracted with Qiagen's Qiacube system and analysed with the Dx CT/NG/MG Assay kit (BioRad).

The concordance of sensitivity and specificity obtained between methods was 98.7 and 100%, respectively. One positive sample gave a negative result for NG with the STD Panel Strip kit. The sample also resulted negative with the NG OligoGen kit (Operon).

Furthermore, the sensitivity of detection of *Chlamydia trachomatis* was analysed with 50 positive samples extracted with the COBAS 4800 platform and analysed with the CT/NG Cobas 4800 kit (Roche).

The concordance of sensitivity obtained between methods was 100%.

Evaluation n° 8

The sensitivity and specificity of detection of *Mycoplasma genitalium* were analysed using 193 DNA samples (86 positive and 107 negative) extracted with Qiagen's Qiacube system and analysed with the Dx CT/NG/MG Assay kit (BioRad).

The concordance of sensitivity and specificity obtained between methods was 94.2 and 100%, respectively.

Evaluation n° 9

The sensitivity of detection of *Ureaplasma spp.* (*Ureaplasma parvum* + *Ureaplasma urealyticum*) was analysed using 31 samples of cervico-vaginal smears preserved in Amies-Stuart medium. These samples had been previously analysed with the Mycoplasma IST 2 kit (BioMerieux), resulting in 30 positives and 1 negative. DNA extraction from the smears was carried out with the DNA Mini Kit (Qiagen).

The concordance of sensitivity obtained between methods was 96.7%. One of the positive samples with the Mycoplasma IST 2 kit proved negative for *Ureaplasma* with the STD Panel Strip kit, but positive for *Mycoplasma hominis*. The analysis of the sample with three individual PCRs with different target regions of the urease gene also gave a negative result for *Ureaplasma*.

The only negative sample from the study resulted positive with the STD Panel Strip kit. It does not apply to assess the concordance of specificity, having studied only one negative sample, but it is important to point out that the Mycoplasma IST 2 kit only gives a positive signal for counts greater than 1,000 CCU/ml, while the STD Panel Strip does not have that limitation. The positive sample was confirmed with individual PCRs of three different regions of the urease gene of *Ureaplasma*.

Furthermore, the sensitivity and specificity of detection of *Ureaplasma spp.* (*Ureaplasma parvum* + *Ureaplasma urealyticum*) were analysed using 64 samples of cervico-vaginal smears preserved in Amies-Stuart medium. These samples had previously been analysed with the Mycoview kit (Zeacon Diagnostics), resulting in 49 positives and 15 negatives. DNA extraction from the smears was carried out with the DNA Mini Kit (Qiagen).

The concordance of sensitivity and specificity obtained between methods was 92 and 73.3%, respectively. Three of the positive samples using the Mycoview kit proved negative for *Ureaplasma* with the STD Panel Strip kit. The analysis of these samples with three individual PCRs with different target regions of the urease gene also gave a negative result for *Ureaplasma* in all cases.

Of the 15 negative samples from the study, four resulted positive with the STD Panel Strip kit. The four positive samples were confirmed with individual PCRs of three different regions of the urease gene of *Ureaplasma*. It should be noted that the Mycoview kit, as with the Mycoplasma IST 2 kit, only gives a positive signal for counts greater than 1,000 CCU/ml, while the STD Panel Strip does not have that limitation.

Evaluation n° 10

The sensitivity and specificity of detection of *Mycoplasma hominis* were analysed using 32 samples of cervico-vaginal smears preserved in Amies-Stuart medium. These samples had been previously analysed with the Mycoplasma IST 2 kit (BioMerieux), resulting in 2 positives and 30 negatives. DNA extraction from the smears was carried out with the DNA Mini Kit (Qiagen).

The concordance of sensitivity and specificity obtained between methods was 50 and 70%, respectively. Of the two positive samples with the Mycoplasma IST 2 kit, one proved negative with the STD Panel Strip, but positive for *Ureaplasma parvum*. This positive was confirmed with individual PCRs of three different regions of the urease gene of *Ureaplasma*. Of the 30 negative samples, 9 resulted positive for *Mycoplasma hominis* with the STD Panel Strip kit. The 9 samples were confirmed as positive with an alternative individual PCR for *Mycoplasma hominis* and hybridisation to a specific probe. As is the case with *Ureaplasma* spp., the Mycoplasma IST 2 kit only gives a positive signal for *Mycoplasma hominis* for counts greater than 1,000 CCU/mL.

Furthermore, the sensitivity and specificity of detection of *Mycoplasma hominis* were analysed using 64 samples of cervico-vaginal smears preserved in Amies-Stuart medium. These samples had previously been analysed with the Mycoview kit (Zeacon Diagnostics), resulting in 19 positives and 45 negatives. DNA extraction from the smears was carried out with the DNA Mini Kit (Qiagen).

The concordance of sensitivity and specificity obtained between methods was 57.9 and 97.8%, respectively. Of the 19 positive samples with the Mycoview kit, 8 proved negative with the STD Panel Strip; all of them were positive for *Ureaplasma*. These negatives were confirmed using an alternative PCR for *M. hominis*; at the same time, the positivity for *Ureaplasma* was confirmed using individual PCRs for three different regions of the urease gene for *Ureaplasma* spp. Of the 45 negative samples, one resulted positive for *Mycoplasma hominis* with STD Panel Strip; that sample was confirmed with an individual alternative PCR for *Mycoplasma hominis* and hybridisation to specific probes. As is the case with the Mycoplasma IST 2 kit, Mycoview only gives a positive result for *M. hominis* for counts greater than 1,000 CCU/ml.

The concordance of sensitivity for *M. hominis* between Mycoplasma IST2 (50%) or Mycoview (57.9%) is very low. These kits take advantage of the biochemical properties of *Ureaplasmas* and *Mycoplasmas* and are based on the observation of a colour change from yellow to purple when urea (*Ureaplasma* spp.) or arginine (*Mycoplasma hominis*) metabolism occurs.

There are publications that indicate that these tests have problems with false positives, particularly for *M. hominis*. In one of these (Biernat-Sudolska M.), an analysis of a total of 5,715 samples is conducted, of which 240 give a positive result for *M. hominis*, while a total of 176 were confirmed as false positives. In all cases, the presence of *Ureaplasmas* with a titre > 10,000 CFU/ml was demonstrated (as occurs in our study).

Evaluation n° 11

The diagnostic sensibility and specificity for the different microorganisms detected with TD Panel Strip was compared with the kits CLART Panel STIsA/STIsB (Genomica). Both kits were used to evaluate 81 DNA samples isolated from smears of different origin (endocervical, vaginal, urethral, pharyngeal, rectal) using the DNA Mini Kit (Qiagen). The obtained results are shown in the following table:

| Final concordance by microorganism between methods | | | | | | | | | | | |
|----------------------------------------------------|------|------|-----|-----|------|------|------|------|------|------|-----|
| | CT | NG | MG | TV | UP | UU | MH | HSV1 | HSV2 | LGV | TP |
| Sensibility | 100 | 94,1 | 100 | 100 | 100 | 100 | 100 | 100 | 100 | NA | 100 |
| Specificity | 96,8 | 100 | 100 | 100 | 89,6 | 88,6 | 95,9 | 100 | 98,5 | 98,8 | 100 |

All the discordant samples were confirmed in favour of STD Panel Strip kit.

POTENTIAL ISSUES

1. No bands appear on the strip, including the control band.

- The conjugate and/or developer were not balanced at room temperature.
- The conjugate and/or developer were not added, or were added in too small quantity.

2. Only the control band appears

- PCR amplification has failed (check with an agarose gel).
- The PCR was not added, or was added in too small quantity in the hybridization step.
- Incorrect hybridization or wash 1 buffers temperature (higher than instructed).
- Incorrect incubator temperature (higher than instructed).

3. The strip has a strong background colour following development

- The washing steps were not performed effectively (wrong timing, buffer quantity too small, cold buffers).

4. Non-homogenous strip staining

- Inadequate shaking (please review the plate thermo-shaker programme).
- The strips were not entirely submerged during the incubations.

5. Unexpected results

- Inadequate incubation and/or buffer/reagent temperatures.
- PCR contamination (verify using a negative control).
- Contamination of adjacent channels owing to the transfer of liquid from one channel to another when adding the hybridization buffer.

The ideal temperature for carrying out the hybridisation process is 42 °C. Between 39 and 42 °C, the test has shown similar performance, without the appearance of non-specific reactions; temperature above 43 °C may cause a slight reduction in band intensity and kit sensitivity.

SENSITIVITY

We define **the limit of detection** of the kit as the minimum amount of DNA/RNA that can be detected with the test.

We define **the test limit** as the minimum amount of DNA/RNA that we can always genotype or always detect with the test.

Both limits have been determined for each one of the pathogens detected with the STD Panel kit. For this, 1/10 serial dilutions of DNA were analysed from each one of them, with three different product batches.

The sensitivity results obtained are shown in the following table:

| | Limit of detection | Limit of the test |
|-----------|--------------------|-------------------|
| CT | 1 copy | 3 copies |
| NG | 2 copies | 6 copies |
| MG | 0,00000128 ng | 0,000032 ng |

| | Limit of detection | Limit of the test |
|------|--------------------|-------------------|
| TV | 1 copy | 3 copies |
| UU | 0,000032 ng | 0,00016 ng |
| MH | 0,000032 ng | 0,00016 ng |
| HSV1 | 6 copies | 12 copies |
| HSV2 | 2 copies | 6 copies |
| LGV | 12 copies | 12 copies |
| TP | 3 copies | 6 copies |

SPECIFICITY

Through the analysis of a large number of negative samples we can exclude the presence of cross reactions with the typical microbiota of the urogenital tract. In this group of microorganisms are included: *Lactobacillus*, *Streptococcus*, *Corynebacterium*, *Staphylococcus*, *Klebsiella*, *Proteus*, *Mycoplasma*, *Atopobium*, *Peptococcus*, *Petostreptococcus*, *Clostridium*, *Bifidobacterium*, *Propionibacterium*, *Eubacterium*, *Bacteroides*, *Prevotella*.

On the other hand, the specificity of the kit has been evaluated by analysing the DNAs of microorganisms that have been more frequently associated with diseases of the urogenital tract. In addition, we have included in the study the bacteria *C. psittaci*, due to its genetic homology in *Omp* gene with the one of *C. trachomatis*.

The absence of cross-reaction has been checked for the following pathogens:

| Microorganism | Type | Quantity in PCR | Result |
|---------------------------------------|----------------------------------------|-----------------|-------------------------------|
| <i>Chlamydia trachomatis</i> | 434 Strain, ATCC, LGVII | 66500 copies | <i>C. trachomatis</i> |
| <i>Neisseria gonorrhoeae</i> | Clinical isolate | 72500 copies | <i>Neisseria gonorrhoeae</i> |
| <i>Mycoplasma genitalium</i> | NCTC10195 | 0,5 ng | <i>Mycoplasma genitalium</i> |
| <i>Trichomonas vaginalis</i> | Clinical isolate | 74000 copies | <i>Trichomonas vaginalis</i> |
| <i>Ureaplasma urealyticum</i> | NCTC 010112 | 0,5 ng | <i>Ureaplasma urealyticum</i> |
| <i>Mycoplasma hominis</i> | NCTC 10111 | 0,5 ng | <i>Mycoplasma hominis</i> |
| <i>Herpes simplex 1</i> | MacIntyre Strain (ATCC VR-539) | 90000 copies | <i>Herpes simplex 1</i> |
| <i>Herpes simplex 2</i> | MS Strain (ATCC VR-540) | 80000 copies | <i>Herpes simplex 2</i> |
| <i>Treponema pallidum</i> | Nichols strain | 3500 copies | <i>Treponema pallidum</i> |
| <i>E.coli</i> | JM109Strain | 88 ng | Negativo |
| <i>Haemophilus ducreyi</i> | NC10945 | 90 ng | Negativo |
| <i>Streptococcus agalactiae</i> | Stableforth G19 Strain (CECT 183) | 27400 copies | Negativo |
| <i>Staphylococcus Aureus mecA (-)</i> | Subspecies Aureus (CECT 239) | 15 ng | Negative |
| <i>Candida albicans</i> | NCPI 3152 Strain(UKNCC) | 100 ng | Negative |
| <i>Staphylococcus Aureus mecA (+)</i> | Subspecies Aureus (N° 13661) | 95,6 ng | Negative |
| <i>Neisseria meningitidis</i> | Serotype B. NCTC 1026. | 37200 copies | Negative |
| <i>Gardnerella vaginalis</i> | ATCC 14018 | 32400 copies | Negative |
| <i>Chlamydophila psittaci</i> | Cepa 6 BC | 30000 copias | Negative |
| HPV | 6, 11 | 5 ng | Negative |
| HPV | 16, 18, 31, 39, 51, 53, 56, 58, 59, 82 | Real sample | No cross-reaction |

In the case of *T. pallidum*, the kit doesn't distinguish *T. pallidum pallidum* (causing syphilis) from *T. pallidum pertenuis* and *T. paraluisancuniculi*.

T. paraluisancuniculi is a bacterium that infects rabbits and there are no described zoonoses.

T. pallidum pertenuis does infect humans, primarily those under age 15. It causes a tropical disease called pian, which affects the skin, bones, and cartilage. It is not a sexually transmitted disease. It is transmitted by direct (person-to-person), non-sexual contact with the exudate of the lesions of an infected person.

HOOK EFFECT

Increasing DNA quantities of each of the pathogens detected with the kit have been analysed, without observing inhibition of the signal or the appearance of cross-reactions that may be associated with the Hook Effect.

There is therefore no Hook Effect.

INTRA-ASSAY PRECISION

The intra-assay precision of STD Panel has been studied examining, first, 5 sensitivity curves (dilutions of a control DNA consisting of a mixture of DNA from all the pathogens detected by the kit) and on the other side, five replicates of 7 DNA samples positive for one or several pathogens.

Overall, the kit has shown a good intra-assay precision, obtaining the same genotyping for the samples and a similar sensitivity.

INTER-DAY PRECISION

The inter-day precision of STD Panel Strip has been studied examining, first, a sensitivity curve (dilutions of a control DNA consisting of a mixture of DNA from all the pathogens detected by the kit) and on the other side, 7 real DNA samples. Each sensitivity curve was constructed by the same operator, using the same product lot and over five consecutive days.

Overall, the kit has shown a good inter-day precision, obtaining the same genotyping for the samples and a similar sensitivity.

INTER-LAB PRECISION

The inter-lab precision has been studied examining 7 real DNA samples, by duplicate. Each analysis was performed by a different operator, on the same day and using the same product lot. To vary the most of the conditions of use of the kit, each operator used a thermal cycler and different thermoshaker.

Overall, the kit has shown a good inter-lab precision, obtaining the same genotyping for the samples and a similar sensitivity.

INTER-LOT PRECISION

The inter-lot precision of STD Panel Strip kit has been studied by the analysis using three different product lots of a sensitivity curve (dilutions of a control DNA consisting of a mixture of DNA from all the pathogens detected by the kit) and 52 real DNA samples.

The kit has shown a good inter-lab precision. The results obtained for the 52 DNA samples were identical in all cases, and a similar sensitivity.

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1. Biernat-Sudolska M., "The need to verify of positive *Mycoplasma hominis* results obtained using the Mycoplasma IST2 tests". *Journal of Laboratory Diagnostics* (2013), vol 49, pags 5-8.
2. Domeika M. et al. "Guidelines for the laboratory diagnosis of *Chlamydia trachomatis* infections in East European countries". *JEADV* (2009). 23: 1353-1363.
3. Seth-Smith H.MB. et al. " Co-evolution of genomes and plasmids within *Chlamydia trachomatis* and the emergence in Sweden of a new varaint strain". *BMC Genomics* (2009), vol 10, pags 239-249.
4. Sokolovskiy E. et al. "Guidelines for the laboratory diagnosis fo syphilis in East European countries". *JEADV* (2009). 23: 623-632.
5. Domeika M et al. "Guidelines for the laboratory diagnosis of trichomoniasis in East European countries". *JEADV* (2010). 24: 1125-1134.
6. Shipitsyna E. "Guidelines for the laboratory diagnosis of *Mycoplasma genitalium* infections in East European Countries". *Acta Derm. Venereol.* (2010). 90: 461-467.
7. Savicheva A. et al. "Guidelines for laboratory diagnosis of *Neisseria gonorrhoeae* in East-European countries". *Acta Medica Lituanica* (2007). 14: 65-74.
8. Domeika M. et al. "Guidelines for the laboratory diagnosis of genital herpes in eastern European countries". www.eurosurveillance.org.
9. Xiao L. eta al. "Detection and characterization of human *Ureaplasma* species and serovars by Real-Time PCR". *Journal of Clinical Microbiology* (2010). 48 (8): 2715-2723.

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OPERON, S.A. - Camino del Plano, 19 - E-50410 Cuarte de Huerva - (Zaragoza) – SPAIN