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ROTA-ADENO-NORO

FOR PROFESSIONAL USE ONLY

PLEASE READ THESE INSTRUCTIONS CAREFULLY BEFORE USING THE PRODUCT

INTRODUCTION

•Intended use

A positive control to check that the following tests work and occur correctly:

- Simple-Stick Rotavirus
- Simple-Stick Rota-Adeno
- Simple-Stick Norovirus
- Simple Rota-Adeno-Noro

For professional use only. The user needs to have basic sampling and laboratory technique knowledge.

•General information

The product consists of a swab containing a dried preparation of a mix of purified Rotavirus and Adenovirus viral proteins and recombinant virus-like particles (non-infectious) belonging to Norovirus GI.1 (Norwalk) and GII.4 (MD145) genotypes, all diluted in a stabilised buffered solution.

It is a quality control for the aforementioned products to check that the detection process occurs correctly.

This product is exclusively for *in vitro* diagnostic use.

•The testing population

The ROTA-ADENO-NORO dry positive control is intended for use by professionals who perform diagnostic tests using the aforementioned tests to ensure they work correctly.

The Rotavirus, Adenovirus and Norovirus diagnostic test is intended for use on the general population, as anyone of any age is susceptible to infection by these viruses.

•Incidence among the infected population for whom the product is intended

These three viruses are mainly responsible for viral gastrointestinal infections worldwide.

Children aged under 5 years old are at high risk of Norovirus infection. The mean prevalence of this virus is 16%, although this can rise to 22%, for example in South America¹.

Rotavirus infection is also one of the main causes of severe gastroenteritis, particularly among children aged 2 to 5 years old². Vaccine development and their inclusion in national immunisation programmes is helping reduce its incidence. As such, the percentage of hospitalisations caused by Rotavirus is 38.0% in countries where the vaccine is yet to be implemented, while it is down to 23.0% where already implemented³.

When it comes to Adenovirus, its prevalence is lower than Norovirus and Rotavirus, although it remains high. A study carried out in northern France found a 5% prevalence for Adenovirus among children aged 0 to 28 months old⁴. In other poorer countries, reports state that gastroenteritis was found in 1.5% of children in Brazil, 5.1% in Nigeria, and 13% in Guatemala⁵.

• Characteristics of the viruses

Rotaviruses are double-stranded RNA viruses belonging to the Reoviridae family. These viruses have a low infectious dose, with transmission via direct contact between one person and another, either through faeces or orally, and less commonly, through contaminated food and water. This virus is classified into seven antigenic serogroups (A to G). Only groups A, B and C infect human beings, with group A the source of almost all cases⁶⁻⁷.

Noroviruses are a type of single-stranded, positive-sense RNA virus belonging to the Caliciviridae family⁸⁻¹⁰. They are highly contagious viruses and their main transmission routes are person-person contact and through contaminated food/water. They are grouped into five genogroups (GI to GV) and they are also classified by genotypes within each genogroup. Most clinical cases are caused by strains from genogroups I and II, with genotypes GI.1 and GII.4 the most common^{11,12}.

Adenoviruses are double-stranded DNA viruses with a genome of 35,000 base pairs contained in an icosahedral capsid. To date, over 100 human Adenoviruses have been characterised, classified into seven species (A-G). The F species only contains two Adenoviruses, 40 and 41, which are responsible for gastrointestinal-type symptomatic conditions¹³.

THE BASIS OR BASIC PRINCIPLES OF THE PRODUCT

The product consists of a swab containing a previously dried preparation comprising a mix of purified Rotavirus and Adenovirus viral proteins, and recombinant virus-like particles (non-infectious) belonging to Norovirus GI.1 (Norwalk) and GII.4 (MD145) genotypes, all diluted in a stabilised buffered solution. The antigen is released when the swab is immersed in the sample dilution buffer (SDB) from any of the tests to which this control applies, obtaining a solution ready for use with the immunochromatographic strip.

On applying the solution to the test, the antibodies present on the coloured beads react with the antigen contained in the positive control. This latex beads/antibody/antigen complex migrates along the membrane, and it is also taken up by other antibodies located in the membrane test zone, causing a red line, two red lines, or a red line and a blue line to appear (depending on the relevant test, see the "READING THE RESULTS" section), should the detection occur correctly (see Fig. 1).

THE PRODUCT'S FORMATS

The ROTA-ADENO-NORO dry positive control is available as an individually packaged swab.

These instructions for use apply to any trade reference for the product 9.050.YYY.12.000.

INCLUDED MATERIALS

- A swab containing ROTA-ADENO-NORO (dry positive control).
- Instructions for use.
- A plastic microtube.
- A plastic pipette.

The product's safety sheet contains information about the composition of the reagents. An electronic copy of the sheet can be requested via email to msds@operon.es.

MATERIAL NOT INCLUDED IN THE PRODUCT

This is a complementary product and is an accessory for the products stated in the intended use section. Therefore, it is to be used with the materials they include and in accordance with the indications in the HOW TO USE section.

The following materials may also be required but are not included:

- A stopwatch.

PRECAUTIONS

1. The source of the ROTA-ADENO-NORO dry positive control ensures that no higher risk infectious agents are present (hepatitis, HIV).
2. The original diluent buffer in which the native protein is prepared, contains sodium azide at 0,095% as anti-microbial agent. Avoid direct contact with skin and mucous membranes. Dispose of appropriately.
3. In the same way, such buffer contains material from animal origin (albumins). Every material from animal origin is always a commercial reference, and has the certificate that evidences it is not infectious.
4. Do not store or prepare food, eat, drink, or smoke in the area where the reagents and samples are handled.
5. On completing the task, clean the work surfaces with 10 % bleach and/or 70 % ethanol. Wash hands with soap and water, rubbing them well.
6. Only use the reagents *in vitro*.
7. Do not use any of the kit's components after the expiry dates.
8. If the outer box is damaged, the product can still be used providing that the primary packaging remains intact.
9. Throw the product away if the primary packaging is either broken or impaired.
10. Remove the swab from the aluminium bag when the test is going to be performed to avoid unnecessary excessive exposure of the product to any potentially damaging environmental factors.
11. Please dispose of the used product in accordance with current legislation.
12. Do not discard the kit's outer box until all the content has been used. The box contains essential information about the product's CE marking.
13. Any serious incident involving the product must be reported to the manufacturer and the relevant authority in the member state where the user and/or patient are located.

STORAGE

The kit can be stored at any temperature from +2 to +30 °C.

The expiry date is printed on the aluminium packaging and on the box label.

Once the swab is resuspended, it can be used for as many measurements as the volume permits, as long as performed immediately. Do not store the resulting solution. Using a stored resuspension means misusing the test.

HOW TO USE

1. Using the pipette or dropper included in the kit, add 500 µl of the sample dilution buffer to the microtube supplied with the dry positive control (up to the 0.5 mark).
Should the buffer come packaged in single-use sample vials, break the upper end of a vial (dropper section) and add the drops required to reach the 0.5 mark on the microtube.
2. Remove the dry positive control swab from its packaging.
3. Introduce the swab into the microtube with the buffer and allow to rehydrate, vigorously rotating it inside the liquid for 30 seconds. Remove the swab, pressing it against the tube walls to release all the liquid, before throwing it away.
4. Add the following to the immunochromatographic device using the plastic pipette supplied with the dry positive control:
 - 5 drops of the rehydrated positive control, in the case of the Simple/Stick Rotavirus, Simple/Stick Rota-Adeno, and Simple/Stick Norovirus products, or
 - 7 drops of the rehydrated positive control, in the case of the Simple Rota-Adeno-Noro product.

The sample is to be added drop by drop to the test addition

zone, allowing it to fully absorb between each drop.

5. Without further manipulation, observe the results at 5, 10, or 15 minutes (in accordance with the test indications for which the positive control is used).

READING THE RESULTS

The positive control provides a positive result using the indicated tests.

Simple/Stick Rotavirus

Two lines should appear, one blue and one red.

- A blue line: this is the control line and indicates that the immunochromatographic process developed correctly.
- A red line: this indicates that Rotavirus was detected in the tested sample.

See figure 1.

Simple/Stick Rota-Adeno

Three lines should appear, one green, one blue, and one red.

- A green line: this is the control line and indicates that the immunochromatographic process developed correctly.
- A red line: this indicates that Rotavirus was detected in the tested sample.
- A blue line: this indicates that Adenovirus was detected in the tested sample.

See figure 1.

Simple/Stick Norovirus

Three lines should appear, one blue and two red.

- A blue line: this is the control line and indicates that the immunochromatographic process developed correctly.
- Red lines: this indicates that Norovirus GI (upper red line) and GII (lower red line) were detected in the tested sample.

See figure 1.

Simple Rota-Adeno-Noro

The following combinations of lines should appear:


1. Rota/Adeno strip (left side of the housing): three lines should appear, one green, one red, and one blue.

- A green line: this is the control line and indicates that the immunochromatographic process developed correctly.
- A red line: this indicates that Rotavirus was detected in the tested sample.
- A blue line: this indicates that Adenovirus was detected in the tested sample.

2. Norovirus strip (right side of the housing): three lines should appear, one green and two red.

- A green line: this is the control line and indicates that the immunochromatographic process developed correctly.
- Red lines: this indicates that Norovirus GI (upper red line) and GII (lower red line) were detected in the tested sample.

See figure 1.

Note: the intensity of the resulting lines can differ from the example images.

Any different result to those shown indicates an error while performing the test. Repeat the test, closely following the instructions for use.

LIMITATIONS OF THE PROCEDURE

1. This positive control is designed to check that the test occurs correctly; it does not infer the test's sensitivity.
2. Using the antigens contained in the swab as soon as they are rehydrated in the sample dilution buffer is recommended.
3. Monitoring the reaction time is important. A shorter than indicated reaction time could affect the test result, which would be considered invalid.
4. Adding the correct amount of positive control to the test is very important.

CAPABILITY OF DETECTION

Three batches of positive control were evaluated using the different products for which they can be used. All cases provided a clear positive visual result.

Peak Area values over 200 for Adenovirus and over 400 for Rotavirus and Norovirus were obtained for the products compatible for use with the OPERON strip reader.

THE HOOK OR PROZONE EFFECT

There is evidence that the Simple/Stick Norovirus, Simple/Stick Rotavirus, and Simple/Stick Rota-Adeno products do not present a Hook effect up to a concentration of 40 µg/ml of Norovirus, 65 µg/ml of Rotavirus, and 130 µg/ml of Adenovirus. In the case of the Simple Rota-Adeno-Noro product, there is evidence that a Hook effect does not appear until a concentration of 392 µg/ml of Rotavirus, 16.25 µg/ml of Adenovirus, and 10 µg/ml of Norovirus. The positive control uses a concentration of antigens below those concentrations, so a Hook effect will not occur.

REPEATABILITY

INTRA-ASSAY PRECISION

Two batches of ROTA-ADENO-NORO dry positive control were evaluated. Five swabs from each batch were rehydrated, with the resulting solution analysed in duplicate using the Simple Rota-Adeno-Noro test. The analysis was performed on the same day and by the same operator. The resulting positive signal intensity showed high intra-assay precision for the dry positive control.

REPRODUCIBILITY

INTER-DAY PRECISION

Three swabs from a batch of ROTA-ADENO-NORO dry positive control were rehydrated. The resulting solution was analysed in duplicate using the Simple Rota-Adeno-Noro test. The analysis was repeated over a five day period by the same operator. The results indicated very good inter-day reproducibility for the dry positive control.

INTER-OPERATOR PRECISION

Three operators evaluated two batches of ROTA-ADENO-NORO dry positive control. Each operator rehydrated three swabs from each positive control batch, analysing the resulting solution in duplicate using the Simple Rota-Adeno-Noro test. The results indicated high inter-operator precision for the dry positive control.

INTER-BATCH PRECISION

Three swabs from three different batches of ROTA-ADENO-NORO dry positive control were rehydrated. The resulting solution was analysed in duplicate using the Simple Rota-Adeno-Noro test by the same operator, on the same day. The results indicated high inter-batch reproducibility for the dry positive control.

MODIFICATIONS REGARDING THE PREVIOUS VERSION OF THIS DOCUMENT

The information about the performance parameters of the product is completed. All the information related to the evaluation of the product and the obtained results is detailed in the document.

The structure of the general information provided at the beginning of the document is improved, breaking down the different sections in order to be coherent with the IVDR requirements.

BIBLIOGRAFÍA / BIBLIOGRAPHY

1. Liao Y. et al., Global prevalence of norovirus in cases of acute gastroenteritis from 1997 to 2021: An updated systematic review and meta-analysis. *Microbial Pathogenesis*, 2021, 161,105259.
2. <https://www.who.int/teams/health-product-policy-and-standards/standards-and-specifications/vaccines-quality/rotavirus>.
3. Aliabadi N., et al. Global impact of rotavirus vaccine introduction on rotavirus hospitalisations among children under 5 years of age, 2008–16: findings from the Global Rotavirus Surveillance Network. *The Lancet*, 2019, 7, 893.
4. Tran A. et al. Prevalence of Rotavirus, Adenovirus, Norovirus, and Astrovirus Infections and Coinfections among Hospitalized Children in Northern France. *Journal of Clinical Microbiology*, 2010, 48, 1943.
5. Lee B., et al. Pediatric acute gastroenteritis due to adenovirus 40/41 in low- and middle-income countries. *Curr Opin Infect Dis*. 2020, 33, 398.
6. 1. F. Bon et al. Prevalence of a group A rotavirus, human calicivirus, astrovirus, and adenovirus type 40 and 41 infections among children with acute gastroenteritis in Dijon, France. *J Clin Microbiol*. 1999, Sept, 3055-3058
7. Umesh D. Parashar et al. Global illness and deaths caused by rotavirus disease in children. *Emerging Infectious Diseases*, 2003, vol. 9, No.5, 565-572.
8. Atmar RL and Estes MK. Diagnosis of non-cultivable gastroenteritis viruses, the human caliciviruses. *Clin Microbiol Rev*. 2011, Vol 14, 15-37.
9. Ribes Fernández JM and Buesa Gómez J. Infecciones por

Norovirus. Enferm Infecc Microbiol Clin. 2010, Vol 28, 51-55.

10. Patel MM, Hall AJ, Vinjé J and Parashar UD. Noroviruses: a comprehensive review. J Clin Vir. 2009, Vol 44, 1-8.

11. Buesa J et al. Sequential evolution of Genotype GII.4 Norovirus variants causing gastroenteritis outbreaks from 2001 to 2006 in Eastern Spain. J Med Vir. 2008, Vol 50, 1288-1295.

12. Siebenga J et al. Food-borne viruses in Europe network report: the norovirus GII.4 2006b (for US named Minerva-like, for Japan Kobe034-like, for UK V6) variant now dominant in early seasonal surveillance. Eurosurveillance, Jan-Mar 2008. Vol: 13 (Issues 1-3).

13. Rafie K., et al., The structure of enteric human adenovirus 41— A leading cause of diarrhea in children. Science Advances. 2021, 7, 1.

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OPERON, S.A. - Camino del Plano, 19 -E-50410