CONTROL POSITIVO SECO RSV-FLU A+B-RESPIRADENO

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DRY RSV-FLU A+B-RESPIRADENO POSITIVE CONTROL

INSTRUCTIONS FOR USE

Positive control contains the respiratory syncytial virus detergent and photochemically inactivated, A and B type strains of the Influenza virus inactivated by thimerosal and beta propiolactone treatment and Adenovirus purified viral proteins diluted in a buffer solution.

Positive control has been prepared to be used with the IMMUNOCHROMATOGRAPHIC STRIPS for RSV, Influenza and respiratory Adenovirus detection manufactured by OPERON S.A. The use of the control with any other device or test and the obtained result is EXCLUSIVE RESPONSIBILITY OF THE USER. OPERON S.A. is not responsible of it at all.

Preparation of the Positive Control from the swab:

- Remove the swab from the pouch right at the moment of the analysis.

*NOTE: all the required materials to analyse the dry RSV-Flu A+B-RespirAdeno Positive Control are included in the OPERON RSV, Influenza and respiratory Adenovirus detection kits.

- Add tempered RESP EXTR.01 diluent buffer up to the 0.5 mark of a 1.5 ml eppendorf-type tube.
- Dip the swab into the eppendorf tube and moisten it well with the diluent buffer. Rotate the swab against the walls of the tube for at least 1 minute to release the antigens. Squeeze the impregnated liquid by pressing the swab against the walls of the tube. Remove the swab from the tube and discard it.

For its use with the RSV, Influenza and respiratory Adenovirus strips in Stick format:

- Using the yellow plastic pipette, add 9 drops of the resuspended RSV-Flu A+B-RespirAdeno Positive Control to a well of a flat bottom 96-well microplate. Alternatively, add a larger volume of the positive control to a test tube, so that the stick is in contact with the liquid but the liquid must not be over the arrows.
- Take the stick out of the aluminium pouch and insert it into the well (or test tube) with the arrows pointing downwards.
- Incubate the test for 15 minutes at room temperature. Read the result on the white area.

For its use with the RSV, Influenza and respiratory Adenovirus strips in Simple format:

- Take the reaction device out of the aluminium pouch.
- Add 5 drops of the resuspended RSV-Flu A+B-RespirAdeno Positive Control to the rectangular window of the plastic device marked with an arrow using a yellow plastic pipette.
- Incubate the test for 15 minutes at room temperature. Read the result on central window of the device.

In the strips for the detection of RSV and respiratory Adenovirus, a blue line that corresponds to the test control band (aligned with "C" marked on the cassette) and a red line just below the blue control band (aligned with "T" marked on the cassette) that indicates the presence of RSV or Adenovirus in the sample, must appear on the strip.

In the strip for the detection of Influenza, a blue line that corresponds to the test control band (aligned with "C" marked on the cassette), a red line below the control band (aligned with "T2" marked on the cassette) that indicates the presence of Influenza type B in the sample and a second red line (aligned with "T1" marked on the cassette) that indicates the presence of Influenza type A in the sample, must appear on the strip.

The stability of the positive control is **24 months** if stored between +2°C and +30°C.

Follow always the Good Laboratory Practises when manipulating this preparation. Please read

the safety data sheet of the product.

The origin of the positive control and its inactivation allow to determine that it is not infectious regarding the most risk infectious agents (Hepatitis, HIV). Please read the safety data sheet of the product.

Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking.

DO-0905180 Rev. 1 – 22.06.2020

Expiry date

Manufactured by

Lot number

Read package insert

For in vitro diagnostic use

This product fulfils the requirements of Directive 98/79/EC on in vitro diagnostic medical devices.

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OPERON, S.A. Camino del Plano, 19 E-50410 Cuarte de Huerva (Zaragoza) Spain