CONTROL POSITIVO EN SECO CRYPTO-GIARDIA-ENTAMOEBIA

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The dry positive control was prepared for use with the Stick/Simple Crypto, Stick/Simple Giardia, Stick/Simple Entamoeba, Stick/Simple Crypto-Giardia, and Stick/Simple Crypto-Giardia-Entamoeba products manufactured by OPERON S.A. THE USER IS SOLELY RESPONSIBLE for any use with an alternative device and the subsequent results produced. OPERON S.A. does not accept any responsibility in this event.

Please proceed as follows, depending on the product being used.

Materials included in the kit

The Crypto-Giardia-Entamoeba dry positive control ONLY CONSISTS OF THE SWAB. The user has to provide any other materials indicated in the following explanations.

Procedure for using with the Crypto, Giardia, and Crypto-Giardia tests:

• Stick format:

- Remove the swab from its container immediately before use.
- Add 300 μ l of the dilution buffer at room temperature (included in the Stick Crypto, Stick Giardia, or Stick Crypto-Giardia kits) to a 1.5 ml Eppendorf-type tube.
- Introduce the swab into the same Eppendorf tube and moisten it well with the buffer, rotating the swab against the tube walls for at least 1 minute. With the swab still in the tube, remove it from the liquid and extract as much of the impregnated liquid as possible by pressing it against the tube walls. Dispose of the swab.
- Add the resuspended Crypto-Giardia-Entamoeba positive control to a microplate well or a test tube/vial. If using a 96 well microplate, add \sim 150 μ l per well, but if using a test tube or vial, add the amount required to soak the strip in the liquid, but without it allowing it to pass the arrow marks.
- Remove the strip from its container and introduce it into the well or vial with the arrows pointing downwards. Discard the silica gel.
- Incubate at room temperature for 5 minutes for Stick Crypto and Stick Giardia, or for 10 minutes for Stick Crypto-Giardia. Read the result in the central section of the strip.

• Simple format:

- Remove the swab from its container immediately before use.
- Unscrew the cap of the individual vial (included in the Simple Crypto, Simple Giardia, and Simple Crypto-Giardia kits), and introduce the swab into the sample dilution buffer ensuring it moistens thoroughly. Rotate the swab along the sides of the vial to release the antigen. Squeeze out the impregnated liquid by pressing the swab against the vial walls.
- Remove the swab, close the vial cap tightly, and shake vigorously to ensure a homogeneous mixture.
- Remove the reaction device from the aluminium bag. Discard the desiccant bag as it is only there to prevent humidity from affecting the test and is not used to perform the test.
- Break the upper end of the vial cap.
- Add 4-5 drops to the reaction device sample addition area.

The Crypto and Giardia tests should show a **BLUE** line to indicate the control line and a **RED** line below the control line that represents the *Cryptosporidium parvum* or *Giardia Lamblia* test line, depending on the product being used.

The Crypto-Giardia test should show a **GREEN** line to indicate the control band, a **RED** line under the control band to indicate the *Giardia Lamblia* test line, and a **BLUE** line under the previous line to indicate the *Cryptosporidium parvum* test line.

The procedure for use with the Entamoeba and Crypto-Giardia-Entamoeba tests:

- Remove the swab from its container immediately before use.
- Add 300 μ l of the dilution buffer at room temperature (included in the Stick/Simple Entamoeba and Stick/Simple Crypto-Giardia-Entamoeba kits) to a 1.5 ml Eppendorf-type tube.
- Introduce the swab into the same Eppendorf tube and moisten it well with the buffer, rotating the swab against the tube walls for at least 1 minute. With the swab still in the tube, remove it from the liquid and extract as much of the impregnated liquid as possible by pressing it against the tube walls. Dispose of the swab.

• Stick format:

- Add the resuspended Crypto-Giardia-Entamoeba positive control to a microplate well or a test tube/vial. If using a 96 well microplate, add \sim 150 μ l per well, but if using a test tube or vial, add the amount required to soak the strip in the liquid, but without it allowing it to pass the arrow marks. Remove the strip from its container and introduce it into the well or vial with the arrows pointing downwards. Discard the silica gel.
- Incubate at room temperature for 5 minutes for Stick Entamoeba or 10 minutes for Stick Crypto-Giardia-Entamoeba. Read the result in the central section of the test.

• Simple format:

- Remove the reaction device from the aluminium bag. Discard the desiccant bag as it is only there to prevent humidity from affecting the test and is not used to perform the test.
 - Add 125 μ l of the resuspended Crypto-Giardia-Entamoeba positive control to the round sample addition zone.
 - Incubate at room temperature for 5 minutes for Simple Entamoeba or 10 minutes for Simple Crypto-Giardia-Entamoeba. Read the result in the central section of the test.
- The Entamoeba test should show a **BLUE** line to indicate the control line and a **RED** line below the control line to indicate the *Entamoeba histolytica* test line.
- The Crypto-Giardia-Entamoeba test should show a **PURPLE** line to indicate the control line, a **GREEN** line below the control line to indicate the *Entamoeba histolytica* test line, a **RED** line below the previous line to indicate the *Giardia Lamblia* test line, and a **BLUE** line below the previous line to indicate the *Cryptosporidium parvum* test line.

This dry positive control is stable for 12 months at a temperature between +2°C and +30°C. ALWAYS follow Good Laboratory Practices when handling the controls.

The Crypto-Giardia-Entamoeba positive control contains an inactivated antigen. Positive control has been evaluated against the most risk infectious agents (Hepatitis, HIV) leading to a negative result, but it should still be handled as potentially infectious material. Please read the product safety data sheet.

 Rev. 01 – 17.07.2020

Expiry date

Manufactured by

Lot number

Read the instructions for use

For in vitro diagnostic use

This product fulfils the requirements of Directive 98/79/EC on in vitro diagnostic medical devices.

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OPERON, S.A. Camino del Plano, 19 E-50410 Cuarte de Huerva (Zaragoza) Spain + 34976 503597