

STICK cGMP

One-step immunochromatographic test for the detection of casein glycomacropeptide (cGMP) in milk.



www.operondx.com

STICK cGMP


One-step immunochromatographic test for the detection of casein glycomacropeptide (cGMP) in milk

The test is very sensitive and it is able to detect levels of GMP that could indicate adulterations of 4% whey cheese in milk. Even if milking, storage, transport and processing conditions are optimal, it could be as low as 2% or 1%.

Highlights

- ✓ Reliable results.
- ✓ High sensitivity.
- ✓ Different coloured bands. No errors in interpretation of results.
- ✓ Easy to use.
- ✓ Results in 5 minutes.
- ✓ Sample of liquid and powdered milk.
- ✓ Long shelf-life: 24 months.
- ✓ Most convenient transportation cost: room temperature.
- ✓ Cost-effective.

Procedure



Results

