



Enfermedades Infecciosas y Microbiología Clínica

www.elsevier.es/eimc



Brief report

Direct identification of pneumococcal serotypes in blood cultures by a PCR-reverse-hybridisation technique



Juan Carlos Sanz^{a,b,*}, Ruth de Luis^c, Susana del Río^c, Susana Gamen^c, Emilia Cercenado^d, Maria Angeles Orellana^e, Jose Yuste^{f,g}

^a Unidad de Microbiología Clínica, Laboratorio Regional de Salud Pública, Dirección General de Salud Pública, Comunidad de Madrid, Centro de Especialidades Médicas Vicente Soldevilla, Madrid, Spain

^b CIBER de Epidemiología y Salud Pública (CIBERESP), Madrid, Spain

^c Operon S.A., Cuarte de Huerva, Zaragoza, Spain

^d Servicio de Microbiología, Hospital General Universitario Gregorio Marañón, Madrid, Spain

^e Servicio de Microbiología, Hospital Universitario 12 de Octubre, Madrid, Spain

^f Laboratorio de Referencia de Neumococo, Centro Nacional de Microbiología, Instituto de Salud Carlos III, Majadahonda, Madrid, Spain

^g CIBER de Enfermedades Respiratorias (CIBERES), Madrid, Spain

ARTICLE INFO

Article history:

Received 19 February 2019

Accepted 29 April 2019

Available online 24 June 2019

Keywords:

Streptococcus pneumoniae

Serotyping

PCR

Blood culture

ABSTRACT

Introduction: The objective of this study was to assess the performance of a technique (S. PneumoStrip test) based on PCR followed by reverse strip hybridisation for the detection of *Streptococcus pneumoniae* serotypes directly in blood culture vials.

Methods: One hundred and ten (110) pairs of isolated strains and their corresponding original blood culture vials were studied in parallel. Pure isolated strains were conventionally serotyped using latex agglutination and the Quellung reaction. The S. PneumoStrip test was carried out directly in the original blood culture samples.

Results: In 102 cases (92.7%), results of the serotype obtained by Quellung coincided with their corresponding original blood cultures typed by S. PneumoStrip.

Conclusions: S. PneumoStrip test is a good alternative technique for direct pneumococcal serotyping in blood culture clinical samples.

© 2019 Elsevier España, S.L.U. and Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica. All rights reserved.

Identificación directa de serotipos de neumococo en frascos de hemocultivo mediante una técnica de PCR seguida de hibridación reversa

RESUMEN

Introducción: El objetivo de este estudio fue evaluar el rendimiento de una técnica de PCR seguida de hibridación reversa (S. PneumoStrip test) para su aplicación directa en muestras de frasco de hemocultivo.

Métodos: Se estudiaron 110 cepas aisladas en sangre y sus correspondientes muestras de frasco de hemocultivo. Las cepas fueron serotipadas convencionalmente mediante aglutinación por látex y reacción de Quellung. El test de S. PneumoStrip se realizó directamente en las muestras originales de hemocultivo.

Resultados: En 102 casos (92,7%) el serotipado por Quellung coincidió con el obtenido directamente en su original frasco de hemocultivo.

Conclusiones: S. PneumoStrip constituye una buena alternativa para el serotipado directo en muestras de frasco de hemocultivo.

© 2019 Elsevier España, S.L.U. y Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica. Todos los derechos reservados.

Palabras clave:

Streptococcus pneumoniae

Serotipado

PCR

Hemocultivo

* Corresponding author.

E-mail address: juan.sanz@salud.madrid.org (J.C. Sanz).

<https://doi.org/10.1016/j.eimc.2019.04.010>

0213-005X/© 2019 Elsevier España, S.L.U. and Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica. All rights reserved.

Introduction

Streptococcus pneumoniae serotyping remains critical in the epidemiological surveillance of invasive pneumococcal disease (IPD). The Quellung reaction is the gold standard for serotyping pneumococci. In recent years, some alternative serotyping procedures that include PCR, real-time PCR, sequencing analysis of specific genes, and whole genome sequencing, have been developed.¹ Some of these techniques are reduced to a limited number of detectable serotypes² or require expensive equipment (sequencing analysis).¹ Recently, a commercial assay “S. PneumoStrip test” (Operon S.A., Zaragoza, Spain) for pneumococcal serotyping in strains isolated in culture, based on PCR followed by reverse strip hybridisation, has been reported as easy, sensitive and specific procedure.³ The targets of the test are 76 serotypes (42 individually and 34 in pairs). The objective of this work was to assess the performance of the “S. PneumoStrip test” for the detection of *S. pneumoniae* serotypes directly in blood culture vials.

Materials and methods

One hundred and ten (110) pairs of isolated strains and their corresponding original blood cultures vials, consecutively collected in three public hospitals within the Community of Madrid between July 2017 and July 2018, were studied in parallel (five corresponded to paediatric cases and the rest were from adult patients). Pure isolated strains were conventionally serotyped using the ImmLex™ Pneumotest Kit (Statens Serum Institut, Copenhagen, Denmark), and the Quellung test using pneumococcal factor anti-sera (Statens Serum Institut, Copenhagen, Denmark) was used to complete the full serotyping of the strains.

The S. PneumoStrip test was carried out directly in blood culture samples. For DNA bacterial extraction, 5 mL of blood culture vials were centrifuged at 1000 revolutions per minute (rpm) for 2 min and the supernatants were subsequently re-centrifuged at 10,000 rpm for 3 min. The pellet was used for DNA isolation and purification by using an automated magnetic particle technology with the Qiagen kit (kit EZ1® DSP Virus; Hilden, Germany).

Amplifications were performed according to the S. PneumoStrip manufacturer's recommendations using an ABI 7500 Fast instrument (Applied Biosystems, Foster City, California, USA). The denatured amplicons were placed in three hybridisation nylon strips that include 76 target serotypes among the three strips. Strip A contains the autolysin (*lytA*) and the capsular (*cpsA*) genes and serotypes 1, 3, 4, 5, 6A, 6B, 6C, 6D, 7F/7A, 9A/9V, 14, 18A, 18B/18C, 18F, 19A, 19F, 23A, 23B and 23F. Strip B includes the pneumolysin (*ply*) gene and serotypes 2, 8, 9N/9L, 10A, 10B, 10F/10C, 11A/11D, 11B, 11C, 11F, 12A/46, 12B/44, 12F, 15B/15C, 17F, 20, 22F/22A, 33F/33A and 37. Finally, strip C again includes the *lytA* gene and serotypes 7B, 7C/40, 15A, 15F, 16F, 19B/19C, 21, 25A/25F, 38, 24A, 24B/24F, 31, 32A/32F, 33B/33D, 33C, 35A, 35C, 35F/47F, 41A and 41.³

Results

The serotypes identified by Quellung from the 110 studied strains belonged to 21 single serotypes (8 [30.0%], 3 [14.5%], 9N [6.4%], 15A [4.5%], 23A [3.6%], 24F [2.7%], 4 [1.8%], 10A [1.8%], 11A [1.8%], 13 [1.8%], 22F [1.8%], 23B [1.8%], 24B [1.8%], 35B [1.8%], 35F [1.8%], 6C [1.8%], 9V [1.8%]), 11B [0.9%], 19A [0.9%], 19F [0.9%], and 31 [0.9%], and 4 pairs of serotypes (12B/12F [9.1%], 33F/33A [2.7%], 16F/16A [1.8%], 28F/28A [0.9%] [in these pairs the sera factor cross-reactions cannot allow us to definitively assign a specific serotype microscopically]). Paired strains and blood culture-detected serotypes are shown in Table 1. In 102

cases (92.7%) the serotype results obtained in the direct strain by Quellung and in the corresponding original blood culture with S. PneumoStrip coincided. In 5 cases (4.5%), the serotype detected in the strains by Quellung was not included within the 76 serotypes of the S. PneumoStrip test and, hence, was not identified directly in blood culture. In 3 cases (2.7%), the S. PneumoStrip results in blood culture were negative (no amplification of *lytA*, *cpsA* and *ply* control genes) despite the Quellung reaction identifying serotypes 15A, 9N and 12B/12F in the strains. These 3 strains isolated in agar plates were re-tested by S. PneumoStrip and the obtained results coincided with the Quellung serotype determinations 15A, 9N and 12F. S. PneumoStrip allowed serotype-level identification for some serotypes that cannot be completely differentiated by the Quellung reaction (serotypes 12B/12F classified by S. PneumoStrip as serotype 12F with simultaneously positive bands for 12B/12F/44 and 12F; and serotypes 16F/16A assigned by S. PneumoStrip as 16F by a specific band). Inversely, the serotype 24F identified by Quellung was not differentiated by S. PneumoStrip (24F/24B). The global sensitivity for capsular pneumococci (detection of *lytA*, *cpsA* and *ply* control genes) directly in blood cultures was 107/110 (97.3%). Specificity was 100% confirming an accurate level of serotype discrimination using S. PneumoStrip, with no misidentification directly in blood cultures.

Discussion

To date, up to 99 serotypes of pneumococcus have been described based on the capsular Quellung reaction with commercial rabbit antisera and molecular approaches.^{4,5} Although the Quellung test still remains the conventional standard procedure for pneumococcal serotyping,^{1,4} it is a labour intensive technique that is time-consuming and involves large panels of antisera. Moreover, an additional limitation of the Quellung test is the fact that it requires pneumococcal isolate growth, and therefore, its sensitivity is greatly affected by the culture method. PCR-based serotyping techniques are faster and less laborious than phenotypic methods.⁴ However, they are limited by the number of targets covered. S. PneumoStrip would permit the identification of 76 serotypes according to their capsular polysaccharide biosynthesis (*cps*) locus sequence. It also would be able to detect unencapsulated strains (positive for *lytA* and *ply* but negative for *cpsA* and any serotype band) as well as multiple-serotype infections.⁶ It would be important in colonisation studies because the presence of multiple pneumococcal serotypes is difficult to recognise with conventional culture and the Quellung reaction methods. Nevertheless, PCR serotyping techniques may not be able to distinguish between closely related serotypes when gene sequences vary by only a few mutations.⁴ In this sense, some serotypes are identified together by the S. PneumoStrip Scheme.⁶ In this study, in some cases where the Quellung showed cross-reactions (pairs 12B/12F and 16F/16A), S. PneumoStrip allowed the precise identification of serotype 12F (simultaneously positive bands for 12B/12F/44 and 12F) and 16F (single positive band for 16F). The 13-valent conjugate vaccine (PCV13) covers serotypes 1, 3, 4, 5, 6A, 6B, 7F, 9V, 14, 18C, 19F, 19A, and 23F. Serotypes 22F and 33F are candidates for the new 15-valent pneumococcal conjugate vaccine (PCV15).^{7,8} The recently patented 20-valent pneumococcal conjugate vaccine (PCV20) adds the five serotypes 8, 10A, 11A, 12F and 15B to the PCV15.⁹ All these serotypes are included in the S. PneumoStrip technique, with some (7F, 9V, 22F, 33F, 11A, 15B, 18C) paired to other serotypes of the same serogroup (7A, 9A, 22A, 33A, 11D, 15C, 18B).

The distribution of serotypes observed in this study may reflect the circulation of the pneumococcal population in our region. In

Table 1

Paired serotypes detected in strains by Quellung and directly in blood culture by S. PneumoStrip Test.

Serotype detected in strains by Quellung	Identification as <i>S. pneumoniae</i> by S. PneumoStrip test	Serotype detected in blood cultures by S. PneumoStrip	N	%	% cumulated	Interpretation
8	+	8	33	30.0		Correctly identified to serotype/serogroup level
3		3	16	14.5		
12B/12F		12F	9	8.2		
9N		9N/9L	6	5.5		
15A		15A	4	3.6		
24B		24F/24B	2	1.8		
23A		23A	4	3.6		
33F/33A		33F/33A	3	2.7		
24F		24F/24B	3	2.7		
4		4	2	1.8		
10A		10A	2	1.8		
11A		11A/11D	2	1.8		
16F/16A		16F	2	1.8		
22F		22F/22A	2	1.8		
23B		23B	2	1.8		
6C		6C	2	1.8		
9V		9A/9V	2	1.8		
35F		35F	2	1.8		
31		31	1	0.9		
11B		11B	1	0.9		
19A		19A	1	0.9		
19F		19F	1	0.9	92.7	
13		Negative	2	1.8		Serotypes not included in S. PneumoStrip
35B		Negative	2	1.8		
28F/28A	+	Negative	1	0.9	4.5	
15A	–	Negative ^a	1	0.9		Not identified by S. PneumoStrip to serotype/serogroup level
9N		Negative ^b	1	0.9		
12B/12F		Negative ^c	1	0.9	2.7	
			110	100	100	

^a The serotype Serotype detected in the isolate strain by S. PneumoStrip was 15A.^b The serotype Serotype detected in the isolate strain by S. PneumoStrip was 9N.^c The serotype Serotype detected in the isolate strain by S. PneumoStrip was 12F.

Madrid, during the period 2013–2015, the most frequent serotypes identified were 8, 3, 19A, 22F, 9N, 11A, 23B and 10A respectively.¹⁰ All these serotypes were detected in this work although some were identified by S. PneumoStrip in pairs (9N/9L, 22F/22A, 11A/11D). The non-conjugate covered serotypes 13, 35B and 28F/28A, that were detected by the Quellung test in isolated strains, but not by S. PneumoStrip, are rare (incidence in Madrid from 2013 to 2015 at <0.01 per 100,000).¹⁰ In the current study, we found three blood culture positive hits (2.7%) resulting as negative cases with S. PneumoStrip. This discrepancy was not due to S. PneumoStrip failures (the results in the strains coincided with the Quellung) and might be related to an amplification failure due to low bacterial load in the blood culture media or may be due to a weak hybridisation step limiting the optical reading of reactive bands. However, the global sensitivity reached for capsular pneumococci detection and for serotype identification was high, providing complete specificity. In conclusion, the S. PneumoStrip test is a good alternative approach for direct pneumococcal serotyping in blood culture clinical samples. It would allow to detect the most frequent serotypes circulating in our region.

Financing

The kits of S. PneumoStrip used in this study were provided Operon S.A., Zaragoza, Spain.

Conflict of interest

Disclosure: MS. Ruth de Luis, Ms Susana del Río and Ms Susana Gamen; OPERON S.A. (Zaragoza, Spain).

References

- Varghese R, Jayaraman R, Veeraraghavan B. Current challenges in the accurate identification of *Streptococcus pneumoniae* and its serogroups/serotypes in the vaccine era. *J Microbiol Methods*. 2017;141:48–54.
- Tarragó D, Fenoll A, Sánchez-Tatay D, Arroyo LA, Muñoz-Almagro C, Esteva C, et al. Identification of pneumococcal serotypes from culture-negative clinical specimens by novel real-time PCR. *Clin Microbiol Infect*. 2008;14:828–34.
- Marimón JM, Morales M, Gamen S, Manrique A, Ercibengoa M, Cilla G. A reverse-hybridization test for the identification of 76 pneumococcal serotypes, 42 individually and 34 in pairs. *J Microbiol Methods*. 2017;143:13–6.
- Jauneikaite E, Tocheva AS, Jefferies JM, Gladstone RA, Faust S, Christodoulides M, et al. Current methods for capsular typing of *Streptococcus pneumoniae*. *J Microbiol Methods*. 2015;113:41–9.
- Kjeldsen C, Slott S, Elverdal PL, Sheppard CL, Kapatai G, Fry NK, et al. Discovery and description of a new serogroup 7 *Streptococcus pneumoniae* serotype, 7D, and structural analysis of 7C and 7D. *Carbohydr Res*. 2018;463:24–31.
- Marimón JM, Ercibengoa M, García-Arenzana JM, Alonso M, Pérez-Trallero E. *Streptococcus pneumoniae* ocular infections, prominent role of unencapsulated isolates in conjunctivitis. *Clin Microbiol Infect*. 2013;19:E298–305.
- Sobanjo-ter Meulen A, Vesikari T, Malacaman EA, Shapiro SA, Dallas MJ, Hoover PA, et al. Safety, tolerability and immunogenicity of 15-valent pneumococcal conjugate vaccine in toddlers previously vaccinated with 7-valent pneumococcal conjugate vaccine. *Pediatr Infect Dis J*. 2015;34:186–94.
- Ermlich SJ, Andrews CP, Folkerth S, Rupp R, Greenberg D, McFetridge RD, et al. Safety and immunogenicity of 15-valent pneumococcal conjugate vaccine in pneumococcal vaccine-naïve adults ≥50 years of age. *Vaccine*. 2018. pii:S0264-410X(18)30348-7.

9. Watson WJ, Jodar Martin-Montalvo LP, Isturiz RE, Reinert RR. Immunogenic compositions comprising conjugated capsular saccharide antigens, kits comprising the same and uses thereof. United States Patent Application Publication. Pub. No.: US 2017/0021006 A1. <https://patentimages.storage.googleapis.com/pdfs/92c5d2c0b1754b21d247/US20170021006A1.pdf>.
10. Latasa Zamalloa P, Sanz Moreno JC, Ordobás Gavín M, Barranco Ordoñez MD, Insúa Marisquerena E, Gil de Miguel Á, et al. Trends of invasive pneumococcal disease and its serotypes in the Autonomous Community of Madrid. *Enferm Infecc Microbiol Clin.* 2018;36:612–20.